



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.010

UPDATED FEBRUARY 11, 2025

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 Resolved Issues2

Version Details 3

 Code Stream3

PUBLICATION HISTORY

Revision	Author	Description
February 11, 2025	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of *eRAD RIS 4.2025.010* resolves a TechLive issue where RGB highlight colors were not consistent, a Payment Control calculation issue on manual overrides and a memory leak issue with GUI controls.



Note that *eRAD RIS 4.2025.009* was an internal test release and was not commercialized.

The *4.2025.009* Resolved Issue #36518 has been merged forward into this *4.2025.010* release.

The *4.2025.009* New feature #36188 will be merged forward into a future release and is not included in this *4.2025.010* release.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Integration	36112	Resolves Integration issue where the TechLive RGB highlight colors did not match values used for TechLive Admin in light and dark modes.
Payment Control	36497	Resolved Payment Control issue where manual override values were being incorrectly divided by two before the insurance calculations, resulting in inaccurate values being displayed.
Thick Client GUI	36518	Resolved Thick Client GUI issue with Telerik controls that was causing memory leaks within the application.

SORTED BY CATEGORY AND REDMINE

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

