



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.040

UPDATED AUGUST 4, 2025

TABLE OF CONTENTS

Summary	2
Release Announcement.....	2
Resolved Issues.....	2
New Features.....	2
Feature Details	3
Portals.....	3
Feature #36605 - Portals - Simplified Procedure Picker.....	3
Version Details	5
Code Stream.....	5

PUBLICATION HISTORY

Revision	Author	Description
August 4, 2025	Michelle Mahabir/ Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2025.040 introduces a Portals enhancement that simplifies the procedure picker control, resolves a TechLive error that occurs when opening a digital form, and merges in the content of ERAD RIS 4.2025.038.1 (which resolves a TechLive error that occurs when launching a second TechLive session.).

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
OS Projects	37398	Resolves a TechLive error that occurs when opening a digital form.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Portals	36605	Portals – Simplified Procedure Picker	This enhancement to Portals updates the procedure picker tree control to a button structure to improve usability and user experience.

SORTED BY CATEGORY AND REDMINE

Refer to the FEATURE DETAILS section below for configuration and usage information.

FEATURE DETAILS

Portals

Feature #36605 - Portals - Simplified Procedure Picker

Summary

This enhancement to Portals updates the procedure picker tree control to a button structure to improve usability and user experience.

Background

The current implementation of the Connect Portal’s procedure picker controls displays options in a hierarchical tree structure. While this layout presents a full view of all previously selected nodes, it can be overwhelming and difficult to navigate for some users—especially those unfamiliar with the underlying procedure hierarchy, or those who prefer a more linear, guided experience.

Feature Description

The Procedure Picker has been updated in both the Patient and Referring Portals to follow a question-and-answer format. This new layout replaces the current "tree" display and instead presents only the current "question" and its associated choices. These choices are presented as buttons, with a guiding sentence or prompt positioned above them to provide context and clarity.

Select an exam

Select the type of exam you want to order.
Note: If this is a medical emergency, call 911. If this exam requires an immediate or STAT appointment, [Please contact us.](#)

Referring Modality __picker

CT	CT (COMPUTERIZED TOMOGRAPHY)	DEXA (BONE DENSITY)	MA
MAMMOGRAPHY	MRI (MAGNETIC RESONANCE IMAGING)	NUCLEAR MEDICINE	PET-CT (POSITRON EMISSION TESTING WITH COMPUTED TOMOGRAPHY)
ULTRASOUND (ALSO KNOWN AS SONOGRAPHY)	XR	X-RAY	BODY
HEAD/NECK			

[ADD ANOTHER EXAM](#)

As users make their selections, only the current step in the process is shown on the screen, keeping the interface clean and focused. Previous choices are hidden from view but can be revisited at any time by stepping back through the selection path.

BACK

Select an exam

Select the type of exam you want to order.

Note: If this is a medical emergency, call 911. If this exam requires an immediate or STAT appointment, Please

BACK

ANGIOGRAPHY
(KNOWN AS CTA)

BODY

NOT SURE TYPE OF
DEXA

AXIAL

APPENDICULAR

Configuration Instructions

Service Team assistance is required to enable this feature.

Portals

Changes to applicationsettings.config

To enable this feature:

- Configure `ProcedurePickerStepWorkflow` to `TRUE`
 - Default is `FALSE`.
 - Can be overridden at the Site Level.

RIS Client

Changes to Organization Lookup Table Settings

- Set `ProcedurePickerStepWorkflow` to `TRUE` for desired Practice associated with site group.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

