



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.044

UPDATED AUGUST 26, 2025

TABLE OF CONTENTS

Summary	2
Release Announcement	2
Resolved Issues	2
New Known Limitations	2
Version Details	3
Code Stream	3

PUBLICATION HISTORY

Revision	Author	Description
August 26, 2025	Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2025.044 resolves a Worklists issue where Auto-Open Worklists were not displaying, and merges in the content of hotfix ERAD RIS 4.2025.042.0.1 (which resolves a Mobile Mammogram link redirection issue and addresses scheduling inconsistencies within the Provider Portal).

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Worklists	#37514	Resolved Worklists issue where Auto-Open Worklists were not displaying.

New Known Limitations

The following new Known Limitations were identified with this release:

- **KNOWN LIMITATION #37515 - APPOINTMENT BOOK NOT DISPLAYING IF CONFIGURED TO AUTO OPEN WITH MORE THAN ONE OTHER WORKLIST.**
 - **Issue:** When Appointment Book is configured to auto open with more than one other worklist, it does not display correctly.
 - **Impact:** When RIS is opened, the appointment book tab opens and displays a correct count of items, but no rows are displayed.
 - **Workaround:** Closing and re-opening the worklist resolves the issue.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

