



## eRAD RIS

# RELEASE ANNOUNCEMENT

### Build 4.2025.019

UPDATED MARCH 26, 2025

#### TABLE OF CONTENTS

**Summary** ..... 2

    Release Announcement..... 2

    New Features..... 2

**Feature Details** ..... 3

    Worklists..... 3

        Feature #36001 – Addition of Short Term Follow Up Worklist..... 3

**Version Details** ..... 5

    Code Stream..... 5

#### PUBLICATION HISTORY

Revision	Author	Description
March 26, 2025	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>▪ Commercial release.</li> </ul>

# SUMMARY

## Release Announcement

*This release of eRAD RIS 4.2025.019 introduces a Short Term Follow Up Worklist for BI-RADS 3 tracking.*

## New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Worklists	36001	Addition of Short Term Follow Up Worklist	This enhancement to Worklists introduces a Short Term Follow Up worklist to support the workflow for BI-RADS 3 tracking.

SORTED BY CATEGORY AND REDMINE

Refer to the FEATURE DETAILS section below for configuration and usage information.

# FEATURE DETAILS

## Worklists

### Feature #36001 – Addition of Short Term Follow Up Worklist

#### Summary


This enhancement to Worklists introduces a Short Term Follow Up Worklist to support the workflow for BI-RADS 3 tracking.

#### Background

Previously, the mammo tracking workflow was being leveraged for additional exams. To increase efficiency a Short Term Follow Up worklist was introduced to allow tracking solely for exams designated as BI-RADS 3.

#### Feature Description

To streamline workflow for BI-RADS 3 exams, a dedicated `Short Term Follow Up` worklist has been added to the `Mammo` dropdown, ensuring efficient tracking and management. Exams designated for further review will be added to the worklist for additional evaluation. Once the necessary follow-up actions are completed and documented, the exam will be removed from the worklist.



The workflow remains consistent with the existing process for other BI-RADS categories.

#### Configuration Instructions

No System Administrator actions are necessary to enable this feature; however, optional configuration is available.

#### RIS Client

##### New Short Term Follow Up Worklist

- Exams marked accordingly will fall to this list for additional follow up.
  - Access will need to be granted to have visibility of this worklist.

##### Changes to AccessString Lookup Table Settings

- Set `WL.Birad3` permissions to `Full` to authorized users.
- Grant access permissions as necessary.

The following related settings were added or updated:

Setting	Default	Purpose
WLBirads3	Value=[None Full], Default=[None]	Controls access to the "Short Term Follow-up worklist" menu item of BI-RADS category 3 items.

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:

