



# eRAD RIS

## RELEASE ANNOUNCEMENT

### Build 4.2025.003

UPDATED JANUARY 10, 2025

#### TABLE OF CONTENTS

**Summary** ..... 2

    Release Announcement..... 2

    Resolved Issues..... 2

**Version Details** ..... 3

    Code Stream..... 3

#### PUBLICATION HISTORY

Revision	Author	Description
January 10, 2025	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>Commercial release.</li> </ul>

# SUMMARY

## Release Announcement

*This release of ERAD RIS 4.2025.003 resolves a UI Plugin CD importing issue and merges the content of ERAD RIS 4.2024.126.1 (resolves PACS error occurring when opening studies for dictation).*

## Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
UI Plugins	36416	Resolves UI Plugin issue that caused CDs to import without a WMI file in the directory, leading to locking within RIS.

SORTED BY CATEGORY AND REDMINE

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:

