



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.036

UPDATED JUNE 20, 2025

TABLE OF CONTENTS

Summary	2
Release Announcement.....	2
New Features.....	2
Resolved Issues.....	2
New Known Limitations.....	2
Feature Details	3
OS Projects.....	3
Feature #37163 – Support Default Language Configuration for KIOSK system....	3
Version Details	4
Code Stream.....	4

PUBLICATION HISTORY

Revision	Author	Description
June 20, 2025	Michelle Mahabir/ Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

This release of *ERAD RIS 4.2025.036* resolves a Patient Search issue, introduces an OS Kiosk language system configuration, and merges in the content *ERAD RIS 4.2025.034.0.1* (which resolves a Provider WF Create Order issue where electronic scripts omitted the selected laterality during order creation.) and *ERAD RIS 4.2025.034.0.2* (which resolves a Provider Images issue loading images during high-usage periods.)

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
OS Project	37163	Support Default Language Configuration for KIOSK system	This enhancement to OS Projects adds a system configuration that ensures the default language presented on the Kiosk aligns with the language displayed in eRAD RIS.

Refer to the FEATURE DETAILS section below for configuration and usage information.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Patient Search	36756	Resolved Patient Search issue that caused an error when searching with both Birth Date and First Name.

New Known Limitations

The following new Known Limitations were identified with this release:

- IDENTIFIED BEHAVIOR (FROM #37020) – MORE THAN ONE KEY IMAGES NOT DISPLAYING AND ANNOTATIONS NOT DISPLAYING IN PACS V8 AND V9 PROVIDER PORTALS.
 - ISSUE:** IN SOME CASES, ANNOTATIONS OR CERTAIN KEY IMAGES MAY NOT APPEAR AS EXPECTED IN THE PROVIDER PORTAL.
 - IMPACT:** OCCASIONALLY, PROVIDERS ARE UNABLE TO SEE KEY IMAGES AND ANNOTATIONS WHEN LOADING IMAGES IN THE PROVIDER PORTAL.
 - WORKAROUND:** NO WORKAROUND HAS BEEN IDENTIFIED

FEATURE DETAILS

OS Projects

Feature #37163 – Support Default Language Configuration for KIOSK system

Summary

This enhancement to OS Projects adds a system configuration that ensures the default language presented on the Kiosk aligns with the language displayed in eRAD RIS.

Background

Previously, the configuration was hardcoded in the Kiosk API only.

Feature Description

A new system configuration, `KIOSKSystemDefaultLanguageCode`, has been introduced to define the default language displayed on the Kiosk to match the language set in the RIS system. The value must be an ISO 639-1 language code (e.g., en, es).

Configuration Instructions

System Administrators must complete the following actions to enable this feature:

RIS Client

Changes to SystemConfig Lookup Table

- Enable this feature by setting `KIOSKSystemDefaultLanguageCode` to desired `language code`.
 - The configured values would need to match the i18n (ISO 639-1) la codes.

The following related settings were added or updated:

Setting	Default	Purpose
KIOSKSystemDefaultLanguageCode	en	(Value=ISO language code as string) – The default ISO 639-1 language code used on the KIOSK Updated by #37163

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

