



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.064

UPDATED DECEMBER 18, 2025

TABLE OF CONTENTS

- Summary 3**
 - Release Announcement3
 - New Features3
- Feature Details 4**
 - Database4
 - Feature #37872 - New RIS Configuration setting to uniquely identify the RIS system for SmartOnboarding4
 - Feature #37875 - SmartOnboarding - Support a Configurable Eligibility Workflow4
 - Digital Forms5
 - Feature #37855 - SmartOnboarding - Digital Form layout update to address iFrame display limitations5
 - SmartOnboarding5
 - Feature #37861 - SmartOnboarding - Configurable exam list5
 - Feature #37871 - SmartOnboarding - Support digital signatures on payment receipts5
- Version Details 7**
 - Code Stream7

PUBLICATION HISTORY

Revision	Author	Description
December 17, 2025	Kevin Brooks	▪ Commercial release.
December 18, 2025	Leigh Dymant	▪ Corrections and clarifications.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2025.064 makes several updates to RIS required for SmartOnboarding connectivity, display and future functionality.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Database	37872	New RIS Configuration setting to uniquely identify the RIS system for SmartOnboarding	This enhancement to Database introduces a new SmartOnboarding system configuration setting to enable a unique system identifier on each RIS system.
Database	37875	SmartOnboarding - Support a Configurable Eligibility Workflow	This enhancement to Database introduces a new field that will be used in the SmartOnboarding eligibility workflow to control how insurers without an eligibility connection will be treated.
Digital Forms	37855	SmartOnboarding - Digital Form layout update to address iFrame display limitations	This enhancement to Digital Forms updates the SmartOnboarding workflow to load the questionnaire into a pop-up window due to iPhone compatibility issues.
SmartOnboarding	37861	SmartOnboarding - Configurable exam list	This enhancement allows the exam list for SmartOnboarding to be modified in each RIS system and not be overwritten upon upgrade.
SmartOnboarding	37871	SmartOnboarding - Support digital signatures on payment receipts	This enhancement to SmartOnboarding introduces a new field that will be used in the SmartOnboarding payment workflow to indicate the patient has approved payment.

SORTED BY CATEGORY AND REDMINE

Refer to the FEATURE DETAILS section below for configuration and usage information.

FEATURE DETAILS

Database

Feature #37872 - New RIS Configuration setting to uniquely identify the RIS system for SmartOnboarding

Summary

This enhancement to Database introduces a new SmartOnboarding system configuration setting to enable a unique system identifier on each RIS system.

Feature Description

With this change, a new `KIOSKRISSystemID` SmartOnboarding system configuration setting has been created to uniquely identify each RIS system when establishing connections to other systems.

Configuration Instructions

System Administrators must complete the following actions to enable this feature:

Changes to SystemConfig Lookup Table

- Configure a unique value for `KIOSKRISSystemID` for your installation.

The following related settings were added or updated:

Setting	Default	Purpose
KIOSKRISSystemID	Value=String, Default=Blank	Unique identifier for SmartOnboarding system for connections to TechLive database. Updated by #37872

Feature #37875 - SmartOnboarding - Support a Configurable Eligibility Workflow

Summary

This enhancement to Database introduces a new field that will be used in the SmartOnboarding eligibility workflow to control how insurers without an eligibility connection will be treated.

Feature Description

With this change, a new `selection_instruction` column has been added to the `l_carrier_type` table that will contain a JSON document.

Configuration Instructions

Service Team assistance is required to enable this feature.

Digital Forms

Feature #37855 - SmartOnboarding - Digital Form layout update to address iFrame display limitations

Summary

This enhancement to Digital Forms updates the SmartOnboarding workflow to load the questionnaire into a pop-up window due to iPhone compatibility issues.

Feature Description

Previously, the SmartOnboarding workflow was failing on iPhone as the questionnaire would not correctly load in an iframe after changes from related FEATURE #37602.

With this change when the screen loads in SmartOnboarding:

- The header and footer will be loaded once again.
- The submit button will always be enabled.
- On clicking `Submit`, the page will reload and display a message (at the top of the page):

`If you are not automatically returned to the check-in process in the next 5 seconds, you can close this screen.`

Configuration Instructions

No System Administrator actions are necessary to enable this feature.

SmartOnboarding

Feature #37861 - SmartOnboarding - Configurable exam list

Summary

This enhancement to Database adds the ability for a site (DB) to override the view for any custom filters to return a custom exam list of active order/studies, necessary for the SmartOnboarding workflow.

Feature Description

With this change, an update to the `c_Kiosk_ExamItemList` stored procedure, will perform an override check to run any custom `z_vw_Kiosk_ExamItemList` view if it exists. Otherwise, the default `ris_connect.c_vw_Kiosk_ExamItemList` will be used.

Configuration Instructions

Service Team assistance is required to enable this feature.

Feature #37871 - SmartOnboarding - Support digital signatures on payment receipts

Summary

This enhancement to SmartOnboarding introduces a new field that will be used in the SmartOnboarding payment workflow to indicate the patient has approved payment.

Feature Description

With this change, a new `Patient_approved_payment_flag` column has been added to the `c_payment` table to identify when patients have indicated they approve the payment during SmartOnboarding.

This will be indicated in SmartOnboarding as a checkbox on the payment screen, and will then be used to add verbiage to the universal receipt or a similar documents.

There is no change to the RIS UI.

Configuration Instructions

No System Administrator actions are necessary to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

