



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2026.001

UPDATED JANUARY 13, 2026

TABLE OF CONTENTS

Summary 2

- Release Announcement2
- New Features2
- Resolved Issues2

Feature Details 2

- Database3
 - Feature #32793 - Validate SQL Server 2022.....3
- Lookup Tables3
 - Feature #37927 - Billing Code Not Updating with Vial Size.....3

Version Details 4

- Code Stream4

PUBLICATION HISTORY

Revision	Author	Description
January 13, 2026	Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

This release of *ERAD RIS 4.2026.001* addresses multiple issues related to recent changes to handling of contrast amounts with decimal values, validates RIS support for Microsoft SQL Server 2022, and merges in the content of hotfix release *ERAD RIS 4.2025.061.1* (which resolves an Insurance Eligibility issue where eligibility details were lost from Worklist Flags and IVT Worklist).

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Database	32793	Validate SQL Server 2022	This enhancement to Database validates RIS support for Microsoft SQL Server 2022.
Lookup Tables	37927	Billing Code Not Updating with Vial Size	This enhancement to Lookup Tables updates the logic for contrast billing calculations to accommodate recent updates to support values that contain decimals.

SORTED BY CATEGORY AND REDMINE

Refer to the [FEATURE DETAILS](#) section below for configuration and usage information.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Billing	34363	Resolved Billing issue where billing codes configured with contrast initially displayed with only one unit in Amount to Collect dialog when scheduled, but updated to vial size or amount used when opened from View/Edit and Perform Exam.
Patient Portal	37933	Resolved Patient Portal issue where sourcecode was being displayed on the "reports" tab and "images" tab when there was no data to display.
View Edit	37938	Resolved View Edit issue where an error dialog appeared when opening a signed study that used multiple contrasts from View/Edit.

SORTED BY CATEGORY AND REDMINE

New Known Limitations

The following new Known Limitations were identified with this release:

- **BUG #37977 - ORDERS SCHEDULED IN PORTAL HAVE NO REASON FOR ORDER WHEN OPENED IN REGISTRATION**
 - **Issue:** Viewing the RIS [Registration](#) screen for an order scheduled in either the Provider or Patient Portal is not displaying the [Reason for Order](#) text.
 - **Impact:** The [Reason for Order](#) on both the [Exam](#) tab and the [Order](#) tab appears blank.
 - **Workaround:** The value that was entered in the portal is still visible on the portal.

FEATURE DETAILS

Database

Feature #32793 - Validate SQL Server 2022

Summary

This enhancement to Database validates RIS support for Microsoft SQL Server 2022.

Feature Description

RIS has been validated to be fully functional when configured with Microsoft SQL Server 2022 database.

No compatibility issues have been identified.

Configuration Instructions

Service Team assistance is required to enable this feature.

Lookup Tables

Feature #37927 - Billing Code Not Updating with Vial Size

Summary

This enhancement to Lookup Tables updates the logic for contrast billing calculations to accommodate recent updates to support values that contain decimals.

Feature Description

Previously Feature #35126 added the ability to enter values that contain decimals into the Contrast lookup table, to properly report the amount used for some drug types for billing. However, when the contrast used was added with exam details, the vial size was not auto-populating into Billing Codes.

With this change, the logic calculating amount to collect will now use a decimal unit that is rounded up.

Note that when `l_billing_code.charge_for_each_unit_flag = 'N'` the ATC will only every show 1 for units because in that case we can only charge for 1 unit.

Configuration Instructions

No System Administrator actions are necessary to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

