



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2026.004

UPDATED FEBRUARY 9, 2026

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 New Features2

 Resolved Issues2

Feature Details 3

 Patient Portal3

 Feature #35146 – Multi-Visit Self-Scheduling via Patient Portal.....3

Version Details 9

 Code Stream9

PUBLICATION HISTORY

Revision	Author	Description
February 9, 2026	Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

This rapid release of

*This release of **ERAD RIS 4.2026.004** introduces a Patient Portal enhancement that provides patients with the ability to self-schedule multiple visits in a single step.*

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Patient Portal	35146	Multi-Visit Self-Scheduling via Patient Portal	This enhancement to Patient Portal will provide patients with the ability to self-schedule multiple visits through the Patient Portal.

Refer to the **FEATURE DETAILS** section below for configuration and usage information.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Patient WF: Portal Pass Icode	37758	Resolves Patient WF: Portal Pass Icode issue where rescheduling a Procedure Plan Exam via Icode returned an error.

FEATURE DETAILS

Patient Portal

Feature #35146 – Multi-Visit Self-Scheduling via Patient Portal

Summary

This enhancement to Patient Portal provides patients with the ability to self-schedule multiple visits through the Patient Portal in a single step.


Background

Previously, the Patient Portal only allowed patients to self-schedule exams if they were associated with the same order. When multiple orders needed to be scheduled, patients were required to complete the scheduling workflow separately for each order. The process was time-consuming and created a less efficient user experience that discouraged patients from self-scheduling.

Feature Description

With these changes, the patient will have the ability to schedule multiple exams within the Patient Portal even if they are not associated with the same order creating a more efficient workflow.

If a patient has multiple ordered exams, a dialog will now display listing all ordered studies when users select **MAKE AN APPOINTMENT**:




If multiple exams are associated with an order, all exams will be auto-selected by default.

Note, this dialog will *not* display if the patient does not have more than one unscheduled ordered exam.

Schedule Something Else


If the patient does not see the exam they want to schedule, they have the option to select the `Schedule Something Else` link that will direct the patient to the procedure picker where they can manually choose the exam.



Note that this option can be disabled via a new `MultiOrderNewScheduleEnabled` configuration setting.

Exams That Cannot Be Scheduled via the Portal


Some exams ordered for the patient cannot be scheduled through the Portal - for example if the requested procedure is not available at the site. In these situations, these exams will be unavailable for selection, and a configurable warning message will be displayed in red text:



THE ORDERED EXAM CHECKBOX WILL IS DISABLED FOR EXAMS THAT CANNOT BE SCHEDULED VIA THE PORTAL.

Insurance Selection with Multiple Exams

When multiple exams are being scheduled, a new `Insurance` tab provides the patient the ability to select more than one billing method if needed:



If toggle is unchecked, the patient will have the ability to enter a different billing method per exam:

for CT Chest W & Abdomen W

Who is paying?

Use my insurance | **I'll pay for this exam** | Workers' Comp & Personal Injury

Please be prepared to pay the full amount for this exam at the time of the appointment. When possible, the estimated amount will be displayed before you confirm the patient's appointment.

NEXT CANCEL

Appointment Summary

PATIENT
Mahmud, Milan
DOB: 1-01-2000(26y) Male
MRN: 1005690706_Retro

INSURANCE [EDIT](#)

Mahmud Carrier
Policy Number:10010010
Group Name:

Self Pay
Mahmud Carrier
Policy Number:10010010
Group Name:

EXAM 1

EXAM DETAILS
MA DIGITAL MAMMO CB DIAG UNI & US BREAST

EXAM 2

EXAM DETAILS
CT Chest W & Abdomen W

EXAM 3

EXAM DETAILS
CT Abdomen W & Pelvis W & CTA Chest

However, if the, Apply this insurance to all orders without specified insurance information. checkbox is checked, the defaulted insurance on the first exam will be applied to all other exams **if they do not have insurance specified.**

BACK

Select Insurance for CT Abdomen W & Pelvis W & CTA Chest W/Wo

Select how this exam will be covered.

Use my insurance | **I'll pay for this exam** | Workers' Comp & Personal Injury

Please be prepared to pay the full amount for this exam at the time of the appointment. When possible, the estimated amount will be displayed before you confirm the patient's appointment.

Would you like to use the same insurance for all other orders that don't have insurance specified?

NEXT CANCEL

Exam Review

The Exam tab first allows users to confirm or adjust procedure selections:

INSURANCE ATTACHMENTS **EXAM** REVIEW SCHEDULE CONFIRM

BACK

Note: If this is a medical emergency, call 911.

What are you looking for?
 Mammography
 This exam is a result of a callback or follow-up for additional imaging from a recent mammogram
Which side?
 Left
MA DIGITAL MAMMO CB DIAG UNI & US BREAST
REMOVE EXAM

You have reached the maximum number of exams that can be scheduled online.

What are you looking for?
 CT
Where does it hurt?
 Chest
 Implant_Type_Location
CT Chest W & Abdomen W
REMOVE EXAM

You have reached the maximum number of exams that can be scheduled online.

What are you looking for?
 CT
 Chest / Abdomen / Pelvis
Where does it hurt?
 Abdomen
 Clisc

Appointment Summary

PATIENT
Mahmud, Milan
DOB: 1-01-2000(26y) Male
MRN: 1005690706_Retro

INSURANCE EDIT


Mahmud Carrier
Policy Number:10010010
Group Name:

Self Pay
Mahmud Carrier
Policy Number:10010010
Group Name:

EXAM 1
EXAM DETAILS
MA DIGITAL MAMMO CB DIAG UNI & US BREAST

EXAM 2
EXAM DETAILS
CT Chest W & Abdomen W

EXAM 3
EXAM DETAILS



The next step on the Exam tab is to review details starting with the Reason for exam:

INSURANCE ATTACHMENTS **EXAM** REVIEW SCHEDULE CONFIRM

BACK

Reason for exam:MA DIGITAL MAMMO CB DIAG UNI & US BREAST

Pain

Reason for exam:CT Chest W & Abdomen W

Pain

Reason for exam:CT Abdomen W & Pelvis W & CTA Chest W/Wo
CT Plan
- XR Facial Bones 1-2 Views
- XR FACIAL BONES MIN 3 VIEWS
- XR Orbits 4 Views
XR SCOLIOSIS STND 2 VW
quicker
- XR RIBS BILAT W PA CHEST
- XR SKULL MIN 4 VIEWS

Pain

For the question **Special Accommodations** answers are selected differently for selected procedures as below
MA DIGITAL MAMMO CB DIAG UNI & US BREAST

Appointment Summary

PATIENT
Mahmud, Milan
DOB: 1-01-2000(26y) Male
MRN: 1005690706_Retro

INSURANCE EDIT


Mahmud Carrier
Policy Number:10010010
Group Name:

Self Pay
Mahmud Carrier
Policy Number:10010010
Group Name:

EXAM 1
EXAM DETAILS
MA DIGITAL MAMMO CB DIAG UNI & US BREAST

EXAM 2
EXAM DETAILS
CT Chest W & Abdomen W

EXAM 3
EXAM DETAILS



Below reason for exam details, pre-screening questions are presented.

When pre-screening answers have been previously entered but differ across orders being scheduled, the patient will be given the opportunity to review the conflicting answers and provide a single response:

For the question **Special Accommodations** answers are selected differently for selected procedures as below

MA DIGITAL MAMMO CB DIAG UNI & US BREAST
 Special Accommodations - No

CT Chest W & Abdomen W:
 Special Accommodations-Yes

CT Abdomen W & Pelvis W & CTA Chest W/Wo:
 Special Accommodations-No

Do you need help with anything?

YES NO

What is your current weight (lbs)?

lbs

Are you diabetic?

YES NO

Have you had any recent tattoos (within the past 14 days), including cosmetics?

YES NO

USERS ARE PROMPTED TO PROVIDE A SINGLE RESPONSE TO SCREENING QUESTIONS.

Completing the Workflow

After the **Exam** tab is completed, the patient proceeds through the workflow's **Review**, **Schedule**, and **Confirm** tabs.

INSURANCE ATTACHMENTS EXAM REVIEW SCHEDULE **CONFIRM**

BACK

Please review the appointment date, time, and location carefully. Click Make an Appointment to schedule your appointment. You can receive an email and a text message confirmation of this appointment. Please remember to bring your prescription with you for your appointment.

EXAM 1 MA DIGITAL MAMMO CB DIAG UNI & US BREAST	7:45 AM	☰
EXAM 2 CT Chest W & Abdomen W	8:05 AM	☰
EXAM 3 CT Abdomen W & Pelvis W & CTA Chest W/Wo	8:40 AM	☰
EXAM 4 CT Plan XR Facial Bones 1-2 Views, XR FACIAL BONES MIN 3 VIEWS, XR Orbits 4 Views	9:10 AM	☰
EXAM 5 XR SCOLIOSIS STND 2 VW	9:25 AM	☰
EXAM 6 quickier XR RIBS BILAT W PA CHEST, XR SKULL MIN 4 VIEWS	9:30 AM	☰

Monday, February 9, 2026
7:45 AM

Fisher
193 Stoner Avenue Suite 200 Westminsters, Maryland 21157
(443) 436-1300

Email Confirmation

Service Team Upgrade Instructions

The Service Team must complete the following actions to deploy this feature:

Portals

See also **CHANGES TO PARAGRAPHCONFIG LOOKUP TABLE.**

Changes to applicationsettings.config

Configure `MultiOrderNewScheduleEnabled` to `TRUE` to display the `Schedule Something Else` link via the pop-up dialog.

- Default is `TRUE`. When `FALSE`, the link will not be displayed on the dialog.

Configuration Instructions

System Administrators must complete the following actions to enable this feature and Service Team assistance is required for some actions:

RIS Client

Changes to ParagraphConfig Lookup Table

The following settings were added:

Setting	Default
PPMultipleOrdersApplyInsuranceForOrdersNotSpecified	Apply this insurance to all orders without insurance information
PPMultipleOTSPopUpAlertMessage	Please uncheck exams that cannot be scheduled.
PPMultipleOTSPopUpCallCenterMessage	This exam cannot be scheduled in the Patient Portal, please call the site to schedule.
PPMultipleOTSPopUpInsuranceAlertMessage	Due to different payment methods, this exam will have to be scheduled separately. Alternatively, please call to schedule.
PPMultipleOTSPopUpNewOrder	Please select the exams you are ready to schedule.
PPMultipleOTSPopUpProviderAlertMessage	Due to different Ordering Providers, this exam will have to be scheduled separately. Alternatively, please call to schedule.
PPMultipleOTSPopUpTitle	Please select the exams you are ready to schedule.

Changes to SystemConfig Lookup Table

- Enable this feature by setting `PPAllowMultipleOrdersToBeSchedule` to `True`


The following related settings were added or updated:

Setting	Default	Purpose
PPAllowMultipleOrdersToBeScheduled	Value=Boolean, Default=[False]	When True, the Patient Portal allows multiple order scheduling. Supports practice level override. Updated by #35146

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:



WARNING: The functionality of some previous releases has not yet been merged forward into this release.

Please carefully review the code stream diagram before deploying to ensure that no functionality is lost when upgrading.