



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2026.007.1

UPDATED MARCH 26, 2026

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 New Features2

Feature Details 3

 Meaningful Use3

 Feature #38201 - CMS Reporting - QRDA III export of CT Quality Measure3

Version Details 4

 Code Stream4

PUBLICATION HISTORY

Revision	Author	Description
March 26, 2026	Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

*This hotfix release of **ERAD RIS 4.2026.007.1** introduces CT Quality Measure report generation as a QRDA III file to support CMS MIPS quality reporting requirements.*

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Meaningful Use	38201	CMS Reporting - QRDA III export of CT Quality Measure	This update to Meaningful Use adds the ability to generate a CT Quality Measure report as a QRDA III file to support CMS MIPS quality reporting requirements.

Refer to the **FEATURE DETAILS** section below for configuration and usage information.

FEATURE DETAILS

Meaningful Use

Feature #38201 - CMS Reporting - QRDA III export of CT Quality Measure

Summary

This update to Meaningful Use adds the ability to generate a CT Quality Measure report as a QRDA III file to support CMS MIPS quality reporting requirements.

Feature Description

With this change, a new CT Quality Measure report can now be generated as a QRDA III file in .xml format to report on all radiologists for a practice.

To generate a report, **Start** and **End** dates, as well as **Practices** and/or **Tax ID** are required.

The screenshot shows a software window titled "Meaningful Use" with a sub-window "Quality Reporting Document (QRDA Category I & III) Export All Patient CCDA". Under the "Measure" section, there are two radio buttons: "Breast Cancer Screening" and "CT Dose Screening". The "CT Dose Screening" radio button is selected and highlighted with a red rectangular box. Below this, the "Criteria" section contains several input fields: "Start" and "End" date pickers (01-01-2025 and 03-31-2025), "Radiologist" and "Practices" dropdown menus, "NPI" and "Tax ID" text boxes (with "12-1234567" entered in the Tax ID field). The "Patient Criteria" section includes "Race", "Insurance" (with "Carrier" selected), "Ethnicity", and "Problems" (with "Problem" selected) dropdown menus. At the bottom of the dialog are buttons for "View / Print", "Export QRDA I", "Export QRDA III", and "Close".

USERS CAN NOW SELECT BETWEEN THE EXISTING BREAST CANCER SCREENING AND NEW CT DOSE SCREENING REPORTS.

Configuration Instructions

No System Administrator actions are necessary to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

