

You can assign a study to a specific Radiologist in eRAD RIS.

From the Perform Exam screen in the Technologist Worklist, in the Associated Radiologists section, click the drop-down for Assigned Radiologist and make a selection.

When the exam is completed in RIS, it will be automatically routed to the Radiologist's Pending Dictation Worklist for the selected Radiologist.

The screenshot shows the eRAD RIS interface for a 'Perform Exam' screen. The window title is 'eRAD RIS (TEST\_NewCity\_TEST as KIM BATHURST [bathurstk])'. The menu bar includes File, Administration, Scheduling, Front Desk, Technologist, Radiologist, Mammo, Transcription, and Help. The current user is NY/HIA. The main area displays a table with columns: Include, Flags, Status, Procedures, Scheduled Date, Referring, Sedation, Room, Accession #, Read STAT, Exam STAT, and Modality Type. A single row is visible: Arrived XR Chest 2 Views Pa and Lateral [71020] - Chest, 10-29-2013 2:06 PM, Refer, Unknown, NKXR1, 1000107462, CR. Below the table are tabs for Patient, MU, Order, Documentation, Billing Codes, Attachments, Notes / Exam Times, and Image Request (0). The 'Documentation' tab is active, showing a 'Tech Notes' section with a text area containing 'This is my notes to the radiologis.' and a 'Copy / Paste' button. To the right of the notes are dropdown menus for 'Primary tech \*' (BATHURST, KIM), 'Assisting tech 1', and 'Assisting tech 2'. Below these are radio buttons for 'Last menstrual period' (mm-dd-yyyy) with options for 'Pregnant' and 'Not pregnant', and checkboxes for 'No longer menstruating', 'Shielded patient', 'Pregnancy test', and 'Blood test'. The 'Procedure' section includes dropdowns for 'Scanner \*' (NKXR1), 'Performed procedure \*' (XR Chest 2 Views Pa and Lateral (71020)), 'Body part \*' (Chest), and 'Change reason code \*'. There is also a 'Change reason note' text area and an 'External report workflow' checkbox. The 'Associated Radiologists' section has three dropdown menus: 'Assigned radiologist', 'Performing radiologist', and 'Assistant radiologist'. The 'Signs & Symptoms' section contains a text area with 'Cough for 2 weeks, chest pain RUQ'. At the bottom are buttons for 'View Images', 'Start', 'Suspend', 'Complete', 'Abort', 'Save', and 'Close'. The status bar at the bottom left shows 'Saved: Walk-in: TEST, Patient #9999999N'.