

The procedure plan lookup table has a “report_together_flag” which determines if studies should have a single linked report. If the procedure plan is set to report together, then one procedure will also need to be set as the primary study.

Procedure Plans:

Procedure Plan Code	Name	Description	Report Together...	Last Updated	Active
MA02	MA DIGITAL DIAGNOSTIC UNILAT MAMMO W BREAST US AND	MA DIGITAL DIAGNOSTIC UNILAT MAMMO W BREAST US AND (composite)	N	11-06-2015 8:24 AM	N
MA03	MA DIGITAL CALLBACK UNILATERAL MAMMO W BREAST US	MA DIGITAL CALLBACK UNILATERAL MAMMO W BREAST US (composite)	N	11-06-2015 8:24 AM	N
MA04	MA DIGITAL CALLBACK BILAT MAMMO W BREAST US	MA DIGITAL CALLBACK BILAT MAMMO W BREAST US (composite)	N	11-06-2015 8:24 AM	N
PPMA02	MA Diagnostic Bilat W Breast US Uni (composite)	MA Diagnostic Bilat W Breast US Ltd Uni (composite)	Y	11-06-2015 8:24 AM	Y
PPMA03	MA Diagnostic Bilat W Breast US Bilat (composite)	MA Diagnostic Bilat W Breast US Ltd Bilat (composite)	Y	11-06-2015 8:24 AM	Y
PPMA04	MA Diagnostic Uni W Breast US Bilat (composite)	MA Diagnostic Uni W Breast US Ltd Bilat (composite)	Y	11-06-2015 8:24 AM	Y
PPMA01	MA Diagnostic Uni W Breast US Uni (composite)	MA Diagnostic Uni W Breast US Ltd Uni (composite)	Y	11-06-2015 8:04 AM	Y

Plan Procedures:

Procedure	Seq	Primary Study Flag	Min Wait Time	Max Wait Time	Duration
G0204 MA Digital Mammo Diag Bilat	0	Y	0	0	20
76642 US Breast Limited, Unilat	1	N	0	0	20

Studies can also be linked in the patient folder.

Study	Procedure	Scheduled Date
Report Dr...	CT STN W Chest W & Abd W [...]	11-22-2012 1:40 PM
Exam Done	CT Soft Tissue Neck Wo [70490...	11-22-2012 1:10 PM
Exam Done	CT Recall for Technical Repeat...	11-22-2012 12:40...
Exam Done	CT Head W/Wo & Sinus W/Wo [...]	11-22-2012 12:10...
Exam Done	CT Chest W/Wo Abdomen W/W...	11-22-2012 11:40...
Exam Done	CT Cervical Spine W/Wo [7212...	11-22-2012 11:10...
Exam Done	CT Abdomen W/Wo & Pelvis W/...	11-22-2012 10:40...
Exam Done	CT Abdomen W & Pelvis W & C...	11-22-2012 10:05...
Exam Done	CT Abdomen W & CTA Chest...	11-22-2012 9:30 AM
Report Dr...	CT 3 Phase [CT13]	11-22-2012 9:00 AM

- Dictate
- View Study
- Schedule From
- IVT
- Edit Patient
- Patient Merge
- View/Edit
- Print Forms
- Attachments
- Portals
- Image Request
- Burn CD
- Update Relevance
- RADAR Quick Message
- View Images
- Report
- Report History
- Send Report
- Distribution History
- Link Studies

Patient: ADAMS, Abigail #1505383DEMO (6) | Pending Dictation WL (42 of 117) | Lookup Tables - AccessString | Reception WL (13) | Orders To Schedule (139) | Schedule C

ags Nuggets Status Procedures Scheduled Date Referring MRN Accession # Confirmation

ags	Nuggets	Status	Procedures	Scheduled Date	Referring	MRN	Accession #	Confirmation
		Arrived	CT Pelvis W [72193] - Pelvis	06-08-2017 4:45 PM	Adams, John	1505383...	1067950	1062679
		Arrived	MA Digital Mammo Diag Bilat w...	06-08-2017 4:20 PM	Adams, John	1505383...	1067948	1062678
		Arrived	US Breast Limited, Bilat [76642...	06-08-2017 4:18 PM	Adams, John	1505383...	1067946	1062677
		Arrived	CT Abdomen W [74160] - Abdo...	06-08-2017 4:15 PM	Adams, John	1505383...	1067951	1062679
		Arrived	MA Digital Mammo Screen Bilat...	06-08-2017 4:00 PM	Adams, John	1505383...	1067949	1062678
		Arrived	MA Digital Mammo Diag Bilat [G...	06-08-2017 3:58 PM	Adams, John	1505383...	1067947	1062677

Link Studies

Linked Studies	Flags	Status	Procedures	Scheduled Date	Locked By	Accession #	Modality T
<input checked="" type="checkbox"/>	\$	Arrived	CT Pelvis W [72193] - Pelvis	06-08-2017 4:45 PM		1067950	CT
<input type="checkbox"/>	\$	Arrived	CT Abdomen W [74160] - Abdomen	06-08-2017 4:15 PM		1067951	CT

Note that linked studies cannot be unlinked.