

HOTFIX

for eRAD RIS
Version 1.0
Build 35
Hotfix 35.005

Table of Contents

1. Purpose..... 3

2. Intended Audience 3

3. Who is affected 3

4. Installing 3

1. Purpose

- 1) Redmine Bug 1152: Mammo exam --> Recommendation codes not saving.
Transcription opens the exam to transcribe but the recommendation field is empty

2. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

3. Who is affected

Installations of eRAD RIS 1.1.35.10140 (RIS application Help About → labeled as 1.1.34.10140)

This is a cumulative HOTFIX. It includes all previous HOTFIXES (35.001 – 35.004)

4. Installing

Please remember to take a backup of all files before replacing with current HOTFIX

To install this hot fix:

- 1) *Replace the contents of the RISDeployment folder with the contents in the ThickClient folder included with this patch. This includes:*
 - a) *Build\YourSystemName\rRIS_1.1.35.10752.zip*
 - b) *Build\rupdate.exe*
 - c) *Core_setup\rRisSetup.msi and setup.exe*
- 2) *Replace the binaries in the rRISService folder with the binaries included with this HOTFIX in the rRISService folder. The will require you to stop the rRISSERVICE*
- 3) *Run the Patch-b35-2.sql file and the rRISCoreData.sql file against your database.*

NOTE other files that are included folders such as MirthConfig were updated in previous HOTFIXES and are also included here. Please refer to previous HOTFIX documentation.

If you have any questions regarding this Hotfix, contact Darcy Noye with the PEI RIS Development Team.