

# Server Update

For eRAD RIS  
Version 2.0  
Build 42  
Update 42.0.8

## **Table of Contents**

1. Defects Resolved .....	3
2. Code Stream.....	4
3. Package Contents .....	5
4. Intended Audience .....	7
5. Who is affected .....	7
6. eRAD RIS Release Version Numbers .....	8
7. Installing .....	9

### **1. Defects Resolved**

Redmine #	Subject	Priority
6196	Tech Status Change is not applied if the front desk has the study open	High

### **Change to Configuration**

A configuration change to the ris.exe.config file was introduced in build 2.42 this setting was designed to lessen the delays when retrieving lookup tables that are larger in size.

This config setting, as mentioned, is located in the ris.exe.config file located in the <appSettings> section of the file

This setting is:

```
<add key="RowLevelCacheDataSetNamesCSV" value="ProcedureCode,TemplateMacro,Indication"/>
```


Where the value can be any lookup table, and a value="\*", will retrieve all lookups.

The setting has resulted in undesired results of possible slowness, depending on the lookup table, and possibly additional bugs. Because of this, we are requesting that all sites comment out this line in the ris.exe.config table. The commented line will look like

```
<!--<add key="RowLevelCacheDataSetNamesCSV" value="ProcedureCode,TemplateMacro,Indication"/> -->
```

## 2. Code Stream

### eRAD RIS - Source Code Evolution



**Legend:**

Light Green = Released software

Gray = Will not be released

Bright Green = Current Release

### **3. Package Contents**

Files included for this server update:

1) RIS\_Service\_Updates contains:

- ABCpdf.dll
- ABCpdf9-32.dll
- ABCpdf9-64.dll
- CDSJobManager.dll
- CDSJobManager.pdb
- CdsLib.dll
- CdsLib.xml
- CdsLib.XmlSerializers.dll
- CommonLib.dll
- CommonLib.pdb
- DDEngine.dll
- DDEngine.pdb
- Interop.ANYMODALCAPTURECTRLLib.dll
- Interop.FAXCOMEXLib.dll
- Microsoft.Practices.EnterpriseLibrary.Caching.dll
- Microsoft.Practices.EnterpriseLibrary.Caching.pdb
- Microsoft.Practices.EnterpriseLibrary.Common.dll
- Microsoft.Practices.EnterpriseLibrary.Common.pdb
- Microsoft.Practices.EnterpriseLibrary.Logging.dll
- Microsoft.Practices.EnterpriseLibrary.Logging.pdb
- Microsoft.Practices.ServiceLocation.dll
- Microsoft.Practices.ServiceLocation.pdb
- Microsoft.Practices.Unity.dll
- Microsoft.Practices.Unity.Interception.dll
- Microsoft.Practices.Unity.Interception.pdb
- Microsoft.Practices.Unity.pdb
- Microsoft.ReportViewer.Common.dll
- Microsoft.ReportViewer.WinForms.dll
- Microsoft.ReportViewer.WinForms.xml
- Newtonsoft.Json.dll
- Newtonsoft.Json.xml
- Ninject.dll
- Ninject.Web.Common.dll
- Ninject.Web.Common.xml
- Ninject.xml
- O2S.Components.PDFRender4NET.dll
- PDFToImages.exe
- rRISServices.exe
- rRISServices.pdb
- rRISServicesx64.exe
- ScheduledTaskManager.dll
- ScheduledTaskManager.pdb
- ServicesConfigWizard.exe
- ServicesConfigWizard.pdb
- SignatureApprovalService.exe
- SignatureApprovalService.pdb
- System.Net.Http.dll
- System.Net.Http.Formatting.dll
- System.Net.Http.Formatting.xml
- System.Net.Http.WebRequest.dll
- System.Net.Http.WebRequest.xml
- System.Net.Http.xml
- System.Web.dll
- System.Web.Helpers.dll
- System.Web.Helpers.xml
- System.Web.Http.dll
- System.Web.Http.SelfHost.dll
- System.Web.Http.SelfHost.xml
- System.Web.Http.WebHost.dll
- System.Web.Http.WebHost.xml
- System.Web.Http.xml
- System.Web.Mvc.dll
- System.Web.Mvc.xml
- System.Web.Optimization.dll
- System.Web.Providers.dll
- System.Web.Razor.dll
- System.Web.Razor.xml
- System.Web.WebPages.Deployment.dll
- System.Web.WebPages.Deployment.xml
- System.Web.WebPages.dll
- System.Web.WebPages.Razor.dll
- System.Web.WebPages.Razor.xml
- System.Web.WebPages.xml
- Turreau.Gilles.FaxDotNet.dll
- Turreau.Gilles.FaxDotNet.pdb
- Turreau.Gilles.FaxDotNet.xml
- WCFService.dll
- WCFService.pdb
- WCFServiceTests.dll
- WCFServiceTests.pdb
- WebApi.dll
- WebApi.pdb
- Wedge.exe
- Wedge.pdb

The rRIS\_Service\_Updates also contains to subfolders:

- A. Deploy folder contains
  - a. Applicationsettings
  - b. ReleaseVersion-connectionstrings
  - c. rRISServices.exe

## **eRAD RIS Server Update**

- d. services
- e. rRISServicesx64.exe

### **B. XSL subfolder contains:**

#### **Subfolder titled “Stylesheets” which contains:**

- a. CDAToHTMLCompleteReport.xslt
- b. CDAToHTMLImpressionSection.xslt
- c. CDAToRIS.xslt
- d. CDAToTxt.xslt
- e. DocumentModelReport.xslt
- f. DocumentModelReportAsXml.xslt

#### **g. CCDA Subfolder contains**

- |                             |                                       |
|-----------------------------|---------------------------------------|
| a. sectionImmunizations.xml | j. header.xml                         |
| b. sectionInstructions.xml  | k. sectionReasonForVisit.xml          |
| c. C-CDA.xml                | l. sectionAdministeredMedications.xml |
| d. sectionMedications.xml   | m. sectionResults.xml                 |
| e. sectionPlanOfCare.xml    | n. sectionAllergies.xml               |
| f. eRadStylesheet.xml       | o. sectionSocialHistory.xml           |
| g. sectionProblemList.xml   | p. sectionEncounters.xml              |
| h. generalTemplates.xml     | q. sectionVitalSigns.xml              |
| i. sectionProcedures.xml    | r. sectionFuncCogStatus.xml           |

### **4. Intended Audience**

This document is created by the RIS Development team for the RadNet/eRAD RIS Implementation teams.

### **5. Who is affected**

Build 42 installs. Server updates are not cumulative. This server update must be applied to 42.0.7

**6. eRAD RIS Release Version Numbers**

Build	Patch	UI Version	Core Version	WS Version	DB Version	Notes
40	17	1.3.40.18780 (3GB)	1.3.40.18756	1.3.40.18756	1.3.40.18646	Core install required for Radiologist User Group. GUI.zip, Core setup, and Web Service updates
40	18	1.3.40.19253(3GB)	1.3.40.18756	1.3.40.18756	1.3.40.19253	GUI.zip, Web Service and database updates
40	19	1.3.40.19909(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	GUI.zip, Web Service
40	20	1.3.40.20086(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip and Recognizer
40	21	1.3.40.20703(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip
40	22	1.3.40.20898(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip
41	-	1.3.41.19135(3GB)	1.3.41.19135	1.3.41.19135	1.3.41.19135	Full version release
41	1	1.3.41.19256(3GB)	1.3.41.19135	1.3.41.19256	1.3.41.19256	GUI.zip, Web Service and database updates
41	2	1.3.41.19620(3GB)	1.3.41.19135	1.3.41.19620	1.3.41.19346	GUI.zip and Web Service
41	2sf	1.3.41.19620SF(3G)	1.3.41.19135	1.3.41.19620	1.3.41.19346	GUI.zip
41	2.1	1.3.41.20622(3GB)	1.3.41.19135	1.3.41.20617	1.3.41.19346	GUI.zip and Web Service
41	3	1.3.41.20125(3GB)	1.3.41.19135	1.3.41.20171	1.3.41.19346	GUI.zip , Web Service and Recognizer
41	4	1.3.41.20285(3GB)	1.3.41.19135	1.3.41.20171	1.3.41.19346	GUI.zip
41	5	1.3.41.20809(3GB)	1.3.41.20708	1.3.41.20792	1.3.41.20738	GUI.zip, Core install, Web Service and database updates
41	6	1.3.41.20912(3GB)	1.3.41.20708	1.3.41.20792	1.3.41.20738	GUI.zip
41.1	-	1.3.41.20918(3GB)	1.3.41.20791	1.3.41.20918	1.3.41.20739	GUI.zip, Core install, Web Service and DB updates, Help files
41.1	1	1.3.41.20987(3GB)	1.3.41.20791	1.3.41.20918	1.3.41.20739	GUI.zip
41.1	2	1.3.41.21162(3GB)	1.3.41.20791	1.3.41.21162	1.3.41.21150	GUI.zip, Web Service and DB updates,
41.1	3	1.3.41.21264(3GB)	1.3.41.20791	1.3.41.21162	1.3.41.21150	GUI.zip
42	-	2.42.0.0.21261(3GB)	2.42.0.0.21261	2.42.0.0.21261	2.42.0.0.21187	Full version release
42	1	2.42.0.1.21597	2.42.0.0.21261	2.42.0.1.21597	2.42.0.0.21187	GUI and Web Service updates
42	2	2.42.0.2.21702	2.42.0.0.21261	2.42.0.2.21702	2.42.0.0.21187	GUI and Web Service updates
42	3	2.42.0.2.21806	2.42.0.0.21261	2.42.0.2.21806	2.42.0.0.21187	GUI and Web Service updates
42	4	2.42.0.2.21862	2.42.0.0.21261	2.42.0.2.21806	2.42.0.4.21870	GUI and Database updates
42	5	2.42.0.5.21924	2.42.0.0.21261	2.42.0.2.21806	2.42.0.4.21870	GUI
42	6	2.42.0.5.22009	2.42.0.0.21261	2.42.0.5.22009	2.42.0.4.21870	GUI and Web Service updates
42	7	2.42.0.7.22090	2.42.0.0.21261	2.42.0.6.22009	2.42.0.4.21870	GUI.zip updates
<b>42</b>	<b>8</b>	<b>2.42.0.7.22090</b>	<b>2.42.0.7.22150</b>	<b>2.42.0.6.22009</b>	<b>2.42.0.4.21870</b>	<b>Web Service updates</b>



### **7. Installing**

- 1) Stop each service listed below. Copy the files from the RIS Service Updates folder to the directory for each rRISService. Be sure to make backup of all files, but especially the rRISServices.exe.config and rRISServices.vshost.exe.config files.

*rRISService directories include:*

- a) rRISService*
- b) CDS rRISService*
- c) Wedge rRISService*
- d) DDE rRISService*

**Restart each service.**

**NOTE:** If your upgrade is coming from 40.016.1 or prior a core install will be required. Versions (builds) of 40.017 and greater do not require a core install as long as you are at least at 40.017 (Core Version 1.3.40.18756)