



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.056

UPDATED JULY 21, 2023

TABLE OF CONTENTS

Summary	2
Release Announcement	2
Resolved Issues	2
New Known Limitations	2
Version Details	3
Code Stream	3

PUBLICATION HISTORY

Revision	Author	Description
July 21, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2023.056 resolves an outbound RADAR communication issue, and addresses a PACS authentication issue.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Provider Images	33660	Addresses a Provider Images issue where users are required to re-authenticate with PACS when their client IP address changes.
RADAR	34170	Resolved RADAR issue where Referring Physician's name was incorrect on emailed reports.

SORTED BY CATEGORY AND REDMINE

New Known Limitations

The following new Known Limitations were identified with this release:

- EXPECTED BEHAVIOR BY DESIGN (FROM #33660) - CLIENT IP CHANGES INVALIDATE PORTAL PACS SESSION
 - **Issue:** Some scenarios remain where Portal users lose PACS authentication when their client IP address changes.
 - **Impact:** In some scenarios, such as switching between Ethernet and Wi-Fi network for communication, users may be prompted to re-authenticate with PACS.
 - **Workaround:** Users may re-authenticate with PACS manually.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

