

# Server Update

For eRAD RIS  
Version 1.0  
Build 41  
Update 41.006

## **Table of Contents**

1. Defects Resolved .....	3
2. Code Stream.....	4
3. Package Contents .....	5
4. Intended Audience .....	5
5. Who is affected .....	5
6. eRAD RIS Release Version Numbers .....	6
7. Installing .....	7

## **RADNET, Inc. – eRAD RIS Server Update**


### **1. Defects Resolved**

<b>Redmine #</b>	<b>Subject</b>	<b>Priority</b>
5584	Maximize Attachments - does not close on deactivation of the reporting screen.	Urgent
5555	Transcription users will get null exceptions when trying to activate ReportTemplateEditor	Urgent
5570	MERGE - Reporting - Transcription user error on first study opened in create report window	High
5581	Spell check on the fly no longer turned on. Squigly lines are not shown under misspelled words	Normal
5559	user is losing practice association and getting error message.	Normal
5388	error on save: Invalid column name 'inverted_face'	Normal

**NOTE:** For bug number 5581 the User Preference “Automatic Spell Check” must be selected for the squiggly lines to appear in the report editor.

## 2. Code Stream

eRad RIS – source code evolution



### **Legend:**

**Green = Released software**

**Red = Will not be released**

**Blue = Future release**

### **3. Package Contents**

Files included for this server update:

- 1) GUI ZIP Updates folder contains
  - a) rRIS\_1.3.41.20912 (3GB).zip

### **4. Intended Audience**

This document is created by the RIS Development team for the RadNet RIS Implementation team.

### **5. Who is affected**

Build 41 installs. This server update must be applied directly to the build 41 core install. Server updates are not cumulative. Any previous server updates must be applied before this one.

## **6. eRAD RIS Release Version Numbers**

<b>Build</b>	<b>Patch</b>	<b>UI Version</b>	<b>Core Version</b>	<b>WS Version</b>	<b>DB Version</b>	<b>Notes</b>
40	10	1.3.40.17435	1.3.40.16659	1.3.40.17435	1.3.40.17370	<b>64 BIT web service.</b> GUI.zip, Webservice.
40	11	1.3.40.17435	1.3.40.16659	1.3.40.17525	1.3.40.17370	Web service update
40	12	1.3.40.17685	1.3.40.16659	1.3.40.17685	1.3.40.17648	GUI.zip, Web service, database updates
40	13	1.3.40.17776	1.3.40.16659	1.3.40.17776	1.3.40.17648	GUI.zip, Web service
40	14	1.3.40.18063(3GB)	1.3.40.16659	1.3.40.18063	1.3.40.18063	GUI.zip, Web service, database updates
40	14.1	1.3.40.18063(3GB)	1.3.40.16659	1.3.40.18107	1.3.40.18063	Web service
40	15	1.3.40.18127(3GB)	1.3.40.16659	1.3.40.18127	1.3.40.18127	GUI.zip, Web service, database updates
40	16	1.1.3.40.18585 (3GB)	1.3.40.16659	1.3.40.18585	1.3.40.18646	GUI.zip, Web service, database updates
40	16.1	1.3.40.18705 (3GB)	1.3.40.16659	1.3.40.18585	1.3.40.18646	GUI.zip
40	17	1.3.40.18780 (3GB)	1.3.40.18756	1.3.40.18756	1.3.40.18646	<b>Core install required</b> for Radiologist User Group. GUI.zip, Core setup, and Web Service updates
40	18	1.3.40.19253(3GB)	1.3.40.18756	1.3.40.18756	1.3.40.19253	GUI.zip, Web Service and database updates
40	19	1.3.40.19909(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	GUI.zip, Web Service
40	20	1.3.40.20086(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip and Recognizer
40	21	1.3.40.20703(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip
40	22		1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip
41	-	1.3.41.19135(3GB)	1.3.41.19135	1.3.41.19135	1.3.41.19135	Full version release
41	1	1.3.41.19356(3GB)	1.3.41.19135	1.3.41.19346	1.3.41.19346	GUI.zip, Web Service and database updates
41	2	1.3.41.19620(3GB)	1.3.41.19135	1.3.41.19620	1.3.41.19346	GUI.zip and Web Service
41	2sf	1.3.41.19620SF(3G)	1.3.41.19135	1.3.41.19620	1.3.41.19346	GUI.zip
41	2.1	1.3.41.20622(3GB)	1.3.41.19135	1.3.41.20617	1.3.41.19346	GUI.zip and Web Service
41	3	1.3.41.20171(3GB)	1.3.41.19135	1.3.41.20171	1.3.41.19346	GUI.zip , Web Service and Recognizer
41	4	1.3.41.20285(3GB)	1.3.41.19135	1.3.41.20171	1.3.41.19346	GUI.zip
41	5	1.3.41.20809(3GB)	1.3.41.20708	1.3.41.20792	1.3.41.20738	GUI.zip, Core install, Web Service and database updates
41	6	1.3.41.20912(3GB)	1.3.41.20708	1.3.41.20792	1.3.41.20738	GUI.zip

## **7. Installing**

- 1) Copy and replace the current rRIS...zip file with the rRIS\_1.3.41.20912.zip file provided with this release.  
Be sure to make a backup of the ris.exe.config file contained within the current .zip folder.