



eRAD RIS

RELEASE ANNOUNCEMENT

Rapid Release 4.2024.008.0.2

UPDATED JANUARY 24, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 New Features2

 Resolved Issues2

Feature Details 3

 Patient WF: Make Appointment3

 Service-Only Feature #35017 – Ability to bypass Insurance step for Mobile Workflow3

Version Details 4

 Code Stream 4

PUBLICATION HISTORY

Revision	Author	Description
January 24, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This rapid release of ERAD RIS 4.2024.008.0.2 addresses a security issue and merges the content of ERAD RIS 4.2024.008.0.1, which introduces the ability to bypass Insurance for Mobile Workflow and resolves a URL site group parameter issue.

New Features

This release introduces the following feature, merged from the uncommercialized ERAD RIS 4.2024.008.0.1 release:

Category	Redmine #	Subject	Description
Patient WF: Make Appointment	35017	Ability to bypass Insurance step for Mobile Workflow	This enhancement to Patient WF: Make Appointment introduces the capability to bypass the insurance step within the Mobile workflow.

Refer to the FEATURE DETAILS section below for configuration and usage information.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Security	35058	This is a service-only change with no customer-facing functionality.

This release also resolves the following issue, merged from the uncommercialized ERAD RIS 4.2024.008.0.1 release:

Category	Redmine #	Subject
Patient WF: Make Appointment	35033	Resolved Patient WF: Make Appointment issue which did not allow the URL Site Group parameter to apply if Site Group was marked inactive.

FEATURE DETAILS

Patient WF: Make Appointment

Service-Only Feature #35017 – Ability to bypass Insurance step for Mobile Workflow

Summary

This enhancement to Patient WF: Make Appointment introduces the capability to bypass the insurance step within the Mobile workflow.

Background

In certain instances, a third party may assume financial responsibility for certain exams so it is not necessary to present the patient with insurance-related inquiries at the time of service.

Previously, the scheduling workflow could be updated to bypass the insurance step via the existing `ScheduleWorkflow` configuration option, but this would disable it for all workflows, not just mobile.

Feature Description

With this change, a new `showInsuranceStep` URL parameter has been created to allow the insurance step to be suppressed within the mobile workflow.

Service Team Upgrade Instructions

The Service Team must complete the following actions to deploy this feature:

Portals

Changes to applicationsettings.config

- To bypass the Insurance step for Mobile Workflow, add the `?showInsuranceStep=False` parameter to the mobile mammo URL.
 - e.g.:

```
http://{patient portal URL}/Mobile/{site group}/{start date}/{end date}?showInsuranceStep=false
```
 - Default is `True` and functionality is unchanged is upon upgrade.
- Alternatively, to configure all workflows (not just mobile workflow) to bypass the insurance step, update the existing `ScheduleWorkflow` configuration to remove `InsuranceStep` from the list of comma-separated values.

Configuration Instructions

Service Team assistance is required to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

