

# **Server Update**

**For eRAD RIS  
Version 1.0  
Build 38  
Update 38.002**

## **Table of Contents**

1. Purpose.....	3
Defects resolved.....	3
2. Package Contents .....	4
3. Intended Audience .....	4
4. Who is affected .....	4
5. eRAD RIS Release Version Numbers .....	4
6. Installing .....	5

## **RADNET, Inc. – eRAD RIS Server Update**

### **1. Purpose**

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update includes the following:

#### **Defects resolved**

Redmine #	Subject	Priority
2498	Close application while minimized, after logging in again the application does not appear	Normal
2246	Reporting - Interactive Dictation SpeechMike buttons not working when set to SubmitAudio	High

### **2. Package Contents**

Files included for this server update:

1. GUI ZIP Updates folder contains
  - a. ris.exe
  - b. ris.pdb

### **3. Intended Audience**

This document is created by the RIS Development team for the RadNet RIS Implementation team.

### **4. Who is affected**

Previous versions of build 38.

Server updates may not be cumulative; therefore this server update must be applied directly to the previous server update.

### **5. eRAD RIS Release Version Numbers**

Build	Patch	UI Version	Core Version	WS Version	DB Version	Notes
37	-	1.2.37.12843	1.2.37.12782	1.2.37.12843	1.2.37.12843	
37	1	1.2.37.13040	1.2.37.13040	1.2.37.13039	1.2.37.13003	
37	2	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	
37	3	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	Just WebAPI updated, no version changes
37	4	1.2.37.13171	1.2.37.13040	1.2.37.13165	1.2.37.13171	
37	5	1.2.37.13347	1.2.37.13040	1.2.37.12979	1.2.37.12979	Bug was discovered with version numbers, this is why it appears we went backwards in the version number but we really did not.
37	6	1.2.37.13665	1.2.37.13040	1.2.37.13580	1.2.37.12979	
37	7	1.2.37.13807	1.2.37.13040	1.2.37.13799	1.2.37.13807	
37	8	1.2.37.13807	1.2.37.13040	1.2.37.13827	1.2.37.13807	
37	9	1.2.37.13807	1.2.37.13040	1.2.37.13839	1.2.37.13807	
37	10	1.2.37.13895	1.2.37.13040	1.2.37.13895	1.2.37.13895	
37	11	1.2.37.14161	1.2.37.13040	1.2.37.14161	1.2.37.14161	
37	12	1.2.37.14285	1.2.37.13040	1.2.37.14300	1.2.37.14285	

## **RADNET, Inc. – eRAD RIS Server Update**

37	13	1.2.37.14326	1.2.37.13040	1.2.37.14326	1.2.37.14326	
38	-	1.2.38.14486	1.2.38.14486	1.2.38.14486	1.2.38.14486	
38	1	1.2.38.14748	1.2.38.14748	1.2.38.14748	1.2.38.14748	
38	2	1.2.38.15182	1.2.38.14748	1.2.38.14748	1.2.38.14748	

## **6. Installing**

- 1) Copy the files in the GUI ZIP Updates folder and replace those in the rRIS\_1.2.38.14748.zip folder. Rename the zip file to be **rRIS\_1.2.38. 15182.zip**

**NOTE:** A core install is not required for this update. Simply exit the application and reopen to get the latest changes.