

Server Update

for eRAD RIS

Version 1.0

Build 36

Update 36.003

Table of Contents

1. Purpose.....	3
2. Intended Audience	4
3. Who is affected	4
4. Installing	4

1. Purpose

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update fixes 2 issues.

- 1.) Redmine Bug 1458 – On Windows XP workstations, the DOB may change for dates that fall between the second Sunday of March to the first Sunday of April and the last Sunday of October to the first Sunday of November. This was due to the change to the start and end of Daylight Savings Time (DST) in 2007. Although there was an XP service pack that allowed XP to properly start and end DST, XP inherently could only support a single algorithm for determining the start and end of DST.
- 2.) Redmine Bug 1505 – Interactive commands not working when using next study workflow.

Files for each reported fix listing about are contained in the 2 folders supplied with this update.

- 1.) Folder GUI ZIP Updates contains the following
 - a. AudioLib.dll
 - b. ris.exe
 - c. WCFClient.exe
- 2.) Folder RIS Service Updates contains the following
 - a. rRISServices.exe
 - b. rRISServices.pdb
 - c. WCFSservice.dll
 - d. WCFSservice.pdb

2. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

3. Who is affected

All previous versions of build 36.

NOTE: Build 36 Server Update 36.002 must be applied before this update.

4. Installing

- 1.) Copy and replace the files **from** the GUI ZIP Updates folder **to** the rRIS_1.1.36.12558.zip file. Rename the .zip to rRIS_1.1.36.12653.zip
- 2.) Copy and replace the files from the RIS Service Updates folder to rRISService directory. Stop each service, replace the files and restart each service.
These directories may include:
 - a. rRISService
 - b. CDS rRISService
 - c. Wedge rRISService
 - d. DDE rRISService