

eRAD RIS

Customer Release Notes

Hotfix Build 3.2018.5.5.6

Updated June 16, 2021

Table of Contents

Summary	3
Intended Audience	3
Who Is Affected	3
Resolved Items	4
Versioning	5
Package Contents	5
Code Stream	5
eRAD RIS Release Version Numbers	6

Publication History

Revision	Author	Description
June 16, 2021	Kevin Brooks / Hilary Saltmarsh	Commercial release.

SUMMARY

Intended Audience

The intended audience for this **CUSTOMER RELEASE NOTES** document is the RIS Administration team for all eRAD RIS customers.



Please carefully review these release notes even if your site will not be upgraded immediately, to identify and communicate any issues that may affect your site.

This document contains information on the purpose, configuration, and operation of new features made available with this release, identifies issues resolved in the release, and highlights any unresolved known limitations.



*There is no separate **SERVICE RELEASE NOTES** edition of this document.*

Who Is Affected

This hotfix release of **eRAD RIS 3.2018.5.5.6** includes two resolved issues, including one that prevented users from logging in to RIS. This hotfix is recommended to be applied on 2018.5.5.X installs.

RESOLVED ITEMS

eRAD RIS

This release resolved the following Redmine tickets:


Redmine #	Subject
30106	Resolved Client Install issue where users were unable to start RIS after an update or when switching between environments.
30168 (Backport of #29530)	Resolved Radiology Reporting issue where RIS could hang when deleting portions of the report.

VERSIONING

Package Contents

	_Documentation	6/15/2021 3:29 PM
	_ReleaseNotes	6/16/2021 4:52 PM
	Client Application	6/15/2021 3:27 PM
	External WebAPI	6/15/2021 3:27 PM
	Identity Service	6/15/2021 3:27 PM
	RIS Service	6/15/2021 3:27 PM
	Service Tools	6/16/2021 8:49 AM
	Build_2018.5.5.6.zip	6/16/2021 8:51 AM
	RISServerMasterCert.pfx	3/31/2016 1:38 PM

Code Stream



Legend:

- Light Green = Previously Released software
- Gray = Internal version, non-release version
- Bright Green = Current Release

eRAD RIS Release Version Numbers

Build	Patch	UI Version	Core Version	WS Version	DB Version	Digital Forms	Patient Portal	UM Portal	Provider Portal	Notes
2018.5	5	3.18.5.5(3GB)	3.18.5.5	3.18.5.5	3.18.5.5.0.02404209	3.18.5.5	3.18.5.5.0.1019	3.18.5.5.0.1019	3.18.5.5.0.1019	Full version release
2018.5	5.1	3.18.5.5.1(3GB)	3.18.5.5.1	3.18.5.5.1	3.18.5.5.1.02460669	3.18.5.5.1	3.18.5.5.1.1023	3.18.5.5.1.1023	3.18.5.5.1.1023	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2018.5	5.2	3.18.5.5.2(3GB)								GUI
2018.5	5.3	3.18.5.5.3(3GB)	3.18.5.5.3	3.18.5.5.3	3.18.5.5.3.02479332					GUI, Web Services, DB
2018.5	5.4	3.18.5.5.4(3GB)	3.18.5.5.4	3.18.5.5.4	3.18.5.5.4.02587388					GUI, Web Services, DB
2018.5	5.5	3.18.5.5.5(3GB)	3.18.5.5.5	3.18.5.5.5	3.18.5.5.5.02615903					GUI, Web Services, DB
2018.5	5.6	3.18.5.5.6(3GB)	3.18.5.5.6	3.18.5.5.6						GUI and Web Services
2018.5	6	3.18.5.6(3GB)	3.18.5.6	3.18.5.6	3.18.5.6.0.02571320	3.18.5.6.0	3.18.5.6.0.1050	3.18.5.6.0.1050	3.18.5.6.0.1050	Full version release
2018.5	6.1	3.18.5.6.1(3GB)	3.18.5.6.1	3.18.5.6.1	3.18.5.6.1.02601339	3.18.5.6.1	3.18.5.6.1.1070	3.18.5.6.1.1070	3.18.5.6.1.1070	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2018.5	6.2	3.18.5.6.2(3GB)	3.18.5.6.2							GUI
2018.5	6.3	3.18.5.6.3(3GB)	3.18.5.6.3	3.18.5.6.3	3.21.5.10.002823110					GUI and Webservices
2020.12.21	-	3.20.12.21.0(3GB)	3.18.5.6	3.20.12.21.0	3.20.12.21.0.02608693	3.20.12.21.0	3.20.12.21.0.897	3.20.12.21.0.897	3.20.12.21.0.897	Full version release
2021.1.4	-	3.21.1.4.0(3GB)	3.21.1.4				3.21.1.4.0.1075	3.21.1.4.0.1075	3.21.1.4.0.1075	GUI and Patient/Provider/UM portals
2021.1.18	-	3.21.1.18.0(3GB)	3.21.1.18	3.21.1.18.0	3.21.1.18.002652234	3.21.1.18.0	3.21.1.18.0.1081	3.21.1.18.0.1081	3.21.1.18.0.1081	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.1	-	3.21.2.1.0(3GB)	3.21.2.1	3.21.2.1.0	3.21.2.1.002672074	3.21.2.1.0	3.21.2.1.0.1090	3.21.2.1.0.1090	3.21.2.1.0.1090	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.15	-	3.21.2.15.0(3GB)	3.21.2.15	3.21.2.15.0	3.21.2.15.002698266	3.21.2.15.0	3.21.2.15.0.1111	3.21.2.15.0.1111	3.21.2.15.0.1111	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.15	1	3.21.2.15.1(3GB)	3.21.2.15.1	3.21.2.15.1		3.21.2.15.1	3.21.2.15.1.1113	3.21.2.15.1.1113	3.21.2.15.1.1113	GUI, Web Services, Patient/Provider/UM portals and Digital Forms
2021.2.15	2	3.21.2.15.2(3GB)	3.21.2.15.2	3.21.2.15.2		3.21.2.15.2	3.21.2.15.2.1122	3.21.2.15.2.1122	3.21.2.15.2.1122	GUI, Web Services, Patient/Provider/UM portals and Digital Forms
2021.2.15	3	3.21.2.15.3(3GB)	3.21.2.15.3							GUI
2021.2.15	4	3.21.2.15.4(3GB)	3.21.2.15.4	3.21.2.15.4	3.21.2.15.402732838	3.21.2.15.4	3.21.2.15.4.1133	3.21.2.15.4.1133	3.21.2.15.4.1133	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.15	5	3.21.2.15.5(3GB)	3.21.2.15.5	3.21.2.15.5	3.21.2.15.502738557					GUI, Web Services and DB
2021.2.15	6	3.21.2.15.6(3GB)	3.21.2.15.6	3.21.2.15.6						GUI and Web Services
2021.3.1	-	3.21.3.1(3GB)	3.21.3.1	3.21.3.1	3.21.3.1.002712308	3.21.3.1	3.21.3.1.0.1119	3.21.3.1.0.1119	3.21.3.1.0.1119	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.3.15	-	3.21.3.15(3GB)	3.21.3.15	3.21.3.15	3.21.3.15.002747230	3.21.3.15	3.21.3.15.0.1144	3.21.3.15.0.1144	3.21.3.15.0.1144	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.4.12	-	3.21.4.12(3GB)	3.21.4.12	3.21.4.12	3.21.4.12.002778929	3.21.4.12	3.21.4.12.0.1167	3.21.4.12.0.1167	3.21.4.12.0.1167	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.5.10	-	3.21.5.10(3GB)	3.21.5.10	3.21.5.10	3.21.5.10.002823110	3.21.5.10	3.21.5.10.0.1189	3.21.5.10.0.1189	3.21.5.10.0.1189	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.5.24	-	3.21.5.24(3GB)	3.21.5.24	3.21.5.24	3.21.5.24.002842120	3.21.5.24	3.21.5.24.0.1205	3.21.5.24.0.1205	3.21.5.24.0.1205	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms