



eRAD RIS

Customer Release Notes

Hotfix Build v3.2021.2.15.5

REVISED MARCH 17, 2021

Table of Contents

Summary	2
Intended Audience	2
Who Is Affected	2
Resolved Issues and Known Limitations	3
Resolved Issues	3
Version Details	4
Package Contents	4
Code Stream	4
eRAD RIS Release Version Numbers	5
Installing	6
Client/GUI	6
Web Service	6
Database	7

Publication History

Revision	Author	Description
March 17, 2021	Kevin Brooks / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Intended Audience

The intended audience for this **Customer Release Notes** document is both the Service and RIS Administration teams for all eRAD RIS customers.

This document describes the purpose, configuration, and operation of new features made available with this release, identifies issues resolved in the release, and highlights any unresolved known limitations.

This information should be used by all sites to determine how these changes affect their site's workflow, and to plan their site's upgrade strategy for eRAD RIS.



*There is no separate **Service Release Notes** edition of this document for this release.*

Who Is Affected

This version contains new features and resolved defects that are recommended to be applied on eRAD RIS VERSION 3.2021.2.15.4.

RESOLVED ISSUES AND KNOWN LIMITATIONS

Resolved Issues















This release resolves the following issues:

Redmine #	Subject
29544	Resolved issue with DB WL cache expiring more often than required. The full fix was not completely included in the 3.2021.2.15.4 hotfix.
29504	Resolved issue with RIS Caching list randomly jumping on worklist refresh.

VERSION DETAILS


Package Contents

The release package includes the following folders:

 _ReleaseNotes	3/17/2021 5:01 PM
 Client Application	3/16/2021 5:54 PM
 DB	3/16/2021 5:54 PM
 External WebAPI	3/16/2021 5:54 PM
 Identity Service	3/16/2021 5:54 PM
 Management Reports	3/16/2021 5:54 PM
 RIS Service	3/16/2021 5:54 PM
 Service Tools	3/16/2021 6:02 PM
 Web Digital Forms	3/16/2021 5:55 PM
 Web Patient Connect	3/16/2021 5:56 PM
 Web Referring Connect	3/16/2021 5:56 PM
 Web UM Connect	3/16/2021 5:56 PM
 Build_2021.2.15.5.zip	3/16/2021 5:57 PM
 RISServerMasterCert.pfx	3/31/2016 1:38 PM

Code Stream

The following source code branches have been merged into this release:



LEGEND:

Light Green = Previously Released software

Gray = Internal version, non-release version

Bright Green = Current Release

eRAD RIS Release Version Numbers

The following table details the version identifiers for components in this release:

Build	Patch	UI Version	Core Version	WS Version	DB Version	Digital Forms	Patient Portal	UM Portal	Provider Portal	Notes
2018.5	1	3.18.5.1.0(3GB)	3.18.5.1.0	3.18.5.1.0	3.18.5.1.0.01916269	3.18.5.1.0	3.18.5.1.0.878	3.18.5.1.0.878	3.18.5.1.0.878	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2018.5	1.1	3.18.5.1.1(3GB)	3.18.5.1.1	3.18.5.1.1	3.18.5.1.1.01983618	3.18.5.1.1	3.18.5.1.1.890	3.18.5.1.1.890	3.18.5.1.1.890	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2018.5	2	3.18.5.2.0(3GB)	3.18.5.2.0	3.18.5.2.0	3.18.5.2.0.02084897	3.18.5.2.0	3.18.5.2.0.924	3.18.5.2.0.924	3.18.5.2.0.924	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2018.5	3	3.18.5.3.0(3GB)	3.18.5.3.0	3.18.5.3.0	3.18.5.3.0.02174049	3.18.5.3.0	3.18.5.3.0.951	3.18.5.3.0.951	3.18.5.3.0.951	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2018.5	4	3.18.5.4.0(3GB)	3.18.5.4.0	3.18.5.4.0	3.18.5.4.0.02214458	3.18.5.4.0	3.18.5.4.0.961	3.18.5.4.0.961	3.18.5.4.0.961	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2018.5	4.1	3.18.5.4.1(3GB)		3.18.5.4.1						GUI and Web Service
2018.5	4.2	3.18.5.4.2(3GB)		3.18.5.4.2						GUI and Web Service
2018.5	4.3			3.18.5.4.3						Web Service
2018.5	4.4	3.18.5.4.4(3GB)	3.18.5.4.4	3.18.5.4.4	3.18.5.4.4.02297855		3.18.5.4.4.982			GUI, Web Service, DB and Patient Portal Updates
2018.5	4.5	3.18.5.4.5(3GB)	3.18.5.4.5	3.18.5.4.5	3.18.5.4.5.02314967	3.18.5.4.5	3.18.5.4.5.988	3.18.5.4.5.988	3.18.5.4.5.988	GUI, Web Service, DB and Patient, Referring and UM Portal Updates
2018.5	4.6	3.18.5.4.6(3GB)	3.18.5.4.6	3.18.5.4.6	3.18.5.4.5.02314967	3.18.5.4.5	3.18.5.4.5.996	3.18.5.4.5.996	3.18.5.4.5.996	GUI, Web Service, Patient, Referring and UM Portal Updates
2018.5	4.7			3.18.5.4.7						Web Services only
2018.5	5	3.18.5.5(3GB)	3.18.5.5	3.18.5.5	3.18.5.5.0.02404209	3.18.5.5	3.18.5.5.0.1019	3.18.5.5.0.1019	3.18.5.5.0.1019	Full version release
2018.5	5.1	3.18.5.5.1(3GB)	3.18.5.5.1	3.18.5.5.1	3.18.5.5.1.02460669	3.18.5.5.1	3.18.5.5.1.1023	3.18.5.5.1.1023	3.18.5.5.1.1023	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2018.5	5.2	3.18.5.5.2(3GB)								GUI
2018.5	5.3	3.18.5.5.3(3GB)	3.18.5.5.3	3.18.5.5.3	3.18.5.5.3.02479332					GUI, Web Services, DB
2018.5	5.4	3.18.5.5.4(3GB)	3.18.5.5.4	3.18.5.5.4	3.18.5.5.4.02587388					GUI, Web Services, DB
2018.5	6	3.18.5.6(3GB)	3.18.5.6	3.18.5.6	3.18.5.6.0.02571320	3.18.5.6.0	3.18.5.6.0.1050	3.18.5.6.0.1050	3.18.5.6.0.1050	Full version release
2018.5	6.1	3.18.5.6.1(3GB)	3.18.5.6.1	3.18.5.6.1	3.18.5.6.1.02601339	3.18.5.6.1	3.18.5.6.1.1070	3.18.5.6.1.1070	3.18.5.6.1.1070	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2020.12.21	-	3.20.12.21.0(3GB)	3.18.5.6	3.20.12.21.0	3.20.12.21.0.02608693	3.20.12.21.0	3.20.12.21.0.897	3.20.12.21.0.897	3.20.12.21.0.897	Full version release
2021.1.4	-	3.21.1.4.0(3GB)	3.21.1.4				3.21.1.4.0.1075	3.21.1.4.0.1075	3.21.1.4.0.1075	GUI and Patient/Provider/UM portals
2021.1.18	-	3.21.1.18.0(3GB)	3.21.1.18	3.21.1.18.0	3.21.1.18.002652234	3.21.1.18.0	3.21.1.18.0.1081	3.21.1.18.0.1081	3.21.1.18.0.1081	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.1	-	3.21.2.1.0(3GB)	3.21.2.1	3.21.2.1.0	3.21.2.1.002672074	3.21.2.1.0	3.21.2.1.0.1090	3.21.2.1.0.1090	3.21.2.1.0.1090	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.15	-	3.21.2.15.0(3GB)	3.21.2.15	3.21.2.15.0	3.21.2.15.002698266	3.21.2.15.0	3.21.2.15.0.1111	3.21.2.15.0.1111	3.21.2.15.0.1111	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.15	1	3.21.2.15.1(3GB)	3.21.2.15.1	3.21.2.15.1		3.21.2.15.1	3.21.2.15.1.1113	3.21.2.15.1.1113	3.21.2.15.1.1113	GUI, Web Services, Patient/Provider/UM portals and Digital Forms
2021.2.15	2	3.21.2.15.2(3GB)	3.21.2.15.2	3.21.2.15.2		3.21.2.15.2	3.21.2.15.2.1122	3.21.2.15.2.1122	3.21.2.15.2.1122	GUI, Web Services, Patient/Provider/UM portals and Digital Forms
2021.2.15	3	3.21.2.15.3(3GB)	3.21.2.15.3							GUI
2021.2.15	4	3.21.2.15.4(3GB)	3.21.2.15.4	3.21.2.15.4	3.21.2.15.402732838	3.21.2.15.4	3.21.2.15.4.1133	3.21.2.15.4.1133	3.21.2.15.4.1133	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.15	5	3.21.2.15.5(3GB)	3.21.2.15.5	3.21.2.15.5	3.21.2.15.502738557					GUI, Web Services and DB
2021.3.1	-	3.21.3.1(3GB)	3.21.3.1	3.21.3.1	3.21.3.1.002712308	3.21.3.1	3.21.3.1.0.1119	3.21.3.1.0.1119	3.21.3.1.0.1119	GUI, Web Services, Patient/Provider/UM portals and Digital Forms

INSTALLING

Client/GUI


Copy and replace the current rRIS...zip file with the eRAD_RIS_2021.2.15.5.zip file provided with this release.

Be sure to make a backup of the ris.exe.config file contained within the current .zip folder.

Web Service

The web service upgrade consists of replacing the files in the RIS Web Service directory with the files supplied in the rRISService folder with this build. Always create a backup of the files being replaced. If using the eRAD_RIS_Backend_Installer.exe, the copying and backup of existing file will be performed for you. This upgrade includes 156 files in the rRISService folder and 17 sub folders. Of those folders, 1 is a folder titled “XSL” that contains a sub folder and files within the subfolder. Another folder is titled “**deploy**” that contains the new .config files that will be deployed with the user of the file called “ServicesConfigWizard.exe” (see instructions below). Folders are (Arabic), en_AU (Australian English), en-ZA (South Africa), he (Hebrew), pt_BR (Brazilian Portuguese), es-AR (Argentina Spanish) and ru-RU (Russian) and are localization folders. The WedgePlugins folder contains various plugin files for RADAR, MIRTH, Ensemble, etc...

Once the files have been copied into place, open each service folder and run the ServicesConfigWizard.exe. All existing configuration settings will be displayed with their current values. If there are any new settings, they will be displayed in **red colored** font with default value applied



Database

To update the database run the database updater tool: RISDataBaseInstaller.exe located in the DB directory of the deployment.