



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.015

UPDATED APRIL 4, 2023

PUBLICATION HISTORY

Revision	Author	Description
April 4, 2023	Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2023.015 resolves various Portal bugs that previously affected daily usage of all portal instances.

New Features

No new features have been introduced with this release.

Resolved Issues

No additional issues were resolved with this release.

This release resolves the following issues:

Category	Redmine #	Subject
Provider Exam Detail Page	27084	Resolved an issue where switching the PACS viewer within the Portal would cause the highlighting border around the thumbnail image to not display.
Admin – User Management	28520	Amended Feature #28520 by adding back in the “Show in Portal Search” checkbox in the Resource Tab.
Provider Search	29929	Resolved a reoccurring error in the SSRS log that was presenting during an advanced exam search within the Provider Portal.
Patient General Display	31137	Resolved a date formatting issue that forced formatting for users that utilize a different language.
Database	32098	Restored the “Referring Portal Activity” core management report as it was returning an error while processing.
Provider Home Screen	32414	Resolved an issue where laterality was missing in exam descriptions within the Provider Portal.
Scheduling	32453	Resolved a row referencing issue within the c_study table that was causing Digital Form answers to not be saved.
Patient Appointment	32620	Resolved an issue that was causing Portal-originated guest appointments to disable the SMS option in the RADAR Quick message pop up in RIS.

New Known Limitations

No new known limitations were identified with this release.