



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.050

UPDATED JULY 6, 2023

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 New Features2

 Resolved Issues2

Feature Details 3

 Insurance Eligibility3

 Feature #33840 - Enhance Group Name criteria to use the co_insurance in the plan vs service type section3

Version Details 5

 Code Stream5

PUBLICATION HISTORY

Revision	Author	Description
July 6, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2023.050 resolves an Insurance Management issue, a Billing issue processing ImaginePay transactions, and introduces new Insurance Eligibility filtering options.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Insurance Eligibility	33840	Group Name criteria to use the co_insurance in the plan vs service type section	This enhancement to Insurance Eligibility updates Eligibility SQL filters to allow exclusion of line items based on source and level values, allowing rules to e.g., exclude all "service_type" level co-insurances when performing Amount to Collect calculations.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Billing	33842	Resolved Billing issue where RIS ImaginePay would timeout when processing Void transactions.
Insurance Management	33722	Resolved Insurance Management issue to handle scenario where individual plan should be used when the family plan shows a zero total deductible for the year.

FEATURE DETAILS

Insurance Eligibility

Feature #33840 - Enhance Group Name criteria to use the co_insurance in the plan vs service type section


Summary

This enhancement to Insurance Eligibility updates Eligibility SQL filters to allow exclusion of line items based on source and level values, allowing rules to e.g., exclude all "service_type" level co-insurances when performing Amount to Collect calculations.

Feature Description

With this change, the Eligibility Special Conditions Table "Edit SQL" filters now allow exclusion of line items based on:

- `Source` - values such as "service_type", "plan"
- `Level` - values such as "Family", "Individual", "Employee"



SQL

```
LineitemType like 'co_ins:%' and Source = 'service_type'
```


Exclude Line Item Flag

	LineItemContents	LineItem Type	Source	Level	Amount	LineItemComment	CarrierCode	MemberId	InsuranceSequenci	Print
▶	{ "percent": "0"...	co_ins: 0	service_type	Individual	0	CAT - PREVENT...	0100demo	#33840	1	010

Line items that will shown

	LineItemContents	LineItem Type	Source	Level	Amount	LineItemComment	CarrierCode	MemberId	InsuranceSequenci	Print
▶	{ "amount": "28...	Deductible : 288...	plan	Family	2883.84	BENEFIT DOES ...	0100demo	#33840	1	C
	{ "amount": "30...	Deductible : 300...	plan	Family	3000.0	BENEFIT DOES ...	0100demo	#33840	1	C
	{ "amount": "48...		plan	Family		AMOUNTS APP...	0100demo	#33840	1	C

This allows flexibility when creating eligibility special conditions selection criteria e.g., support for including all co_insurance that exists at the `plan` level or to exclude all the co_insurance that exists in the `service_type` section, while also supporting combining with other criteria such as the Group Name.



Configuration Instructions

No System Administrator actions are necessary to enable this feature; however, optional configuration is available:

RIS Client

Changes to EligibilitySpecialCondition Lookup Table Settings

- Update rules as desired by modifying the **Matching String** logic via the **Edit SQL...** dialog to add additional **Source** and **Level** criteria.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

