



eRAD RIS

RELEASE ANNOUNCEMENT

Rapid Release 4.2024.017.0.2

UPDATED FEBRUARY 14, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 Resolved Issues2

Version Details 3

 Code Stream3

PUBLICATION HISTORY

Revision	Author	Description
February 14, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> Commercial release.

SUMMARY

Release Announcement

This rapid release of eRAD RIS 4.2024.017.0.2 resolves connection issue to the Epic system.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Provider Images	35037	Resolved Provider Images connection issue by updating Epic decryption process from a computeHash function to a base64ToString call.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

