

Server Update

for eRAD RIS

Version 1.0

Build 37

Update 37.002

Table of Contents

1. Purpose.....	3
2. Intended Audience	3
3. Who is affected	3
4. Installing	4

1. Purpose

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update includes the following:

1. Logged as Redmine bug 1669 - Reporting - Custom Field Exam Title not populating laterality.

Files included are:

- a) WCFSservice.dll
- b) WCFSservice.pdb

2. Logged as Redmine bug 1640 - Can't open image request editing form from Image request WL.

Files included are:

- 1) UILib.dll

2. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

3. Who is affected

All previous versions of build 37. Current build is 1.2.37.13040.

4. Installing

- 1) Stop each service listed below. Copy the files from the RIS Service Updates folder to the directory for each rRISService.


rRISService directories include:

- a) *rRISService*
- b) *CDS rRISService*
- c) *Wedge rRISService*
- d) *DDE rRISService*

Restart each service

- 2) Copy the UILib.dll file from the GUI ZIP Updates folder and replace the file of the same name in the rRIS_1.2.37.13040.zip file.

This may require you to unzip the file, copy and replace the UILib.exe file. Zip the contents again by select all of the files in the rRIS_1.2.37.13040 folder, right click, select Send to > Compressed (zipped) folder.



Rename the .zip file that was created to rRIS_1.2.37.13040a.zip. Move this .zip file up one directory. Remove the old rRIS_1.2.37.13040.zip file and the uncompressed rRIS_1.2.37.13040 folder.

The next time the users log in, updates will automatically be detected because the .zip file has been renamed.