



## eRAD RIS

# RELEASE ANNOUNCEMENT

Build 4.2024.035

UPDATED MARCH 15, 2024

### TABLE OF CONTENTS

**Summary** ..... 2

    Release Announcement .....2

    New Features .....2

    Resolved Issues .....2

**Feature Details** ..... 3

    ICD .....3

    Feature #33388 – Display studies in Reason for Study dialog .....3

**Version Details** ..... 4

    Code Stream .....4

### PUBLICATION HISTORY

Revision	Author	Description
March 15, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>Commercial release.</li> </ul>

# SUMMARY

## Release Announcement

*This release of eRAD RIS 4.2024.035 updates the Reason for Study dialog and resolves a scheduling related issue.*

## New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
ICD	33388	Display studies in Reason for Study dialog	This enhancement to ICD updates the Reason for Study dialog to add a definitive indication of which studies require documentation.

SORTED BY CATEGORY AND REDMINE

Refer to the FEATURE DETAILS section below for configuration and usage information.

## Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Scheduling	35192	Resolved Scheduling issue where an error would occur when attempting to add studies during the scheduling process.

SORTED BY CATEGORY AND REDMINE

# FEATURE DETAILS

## ICD

### Feature #33388 – Display studies in Reason for Study dialog

#### Summary

This enhancement to ICD updates the Reason for Study dialog to add a definitive indication of which studies require documentation.


#### Background

Previously, the dialog alerting that a reason for study was required did not indicate specifically which studies required attention. This was cumbersome to the end users especially if the patient was scheduling multiple studies.

#### Feature Description

The existing dialog has been enhanced to provide the end user with a clear understanding of which studies require a reason for study to be entered.

The dialog will present while saving an order, while scheduling, and during the registration process when necessary.



UPDATED REASON FOR STUDY DIALOG

#### Configuration Instructions

No System Administrator actions are necessary to enable this feature.

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:

