



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.070

UPDATED JUNE 19, 2024

TABLE OF CONTENTS

Summary	2
Release Announcement.....	2
Resolved Issues.....	2
Version Details	3
Code Stream.....	3

PUBLICATION HISTORY

Revision	Author	Description
June 19, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2024.070 resolves an AutoCert error in the View/Edit screen, addresses an error when accessing scanned documents, improves performance within the IVT workflow, and merges in the content of ERAD RIS 4.2024.062.3 (which resolves two issues pertaining to the generation of the C-CDA).

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
IVT/Precert	35591	Resolved IVT/Precert issue where an error would occur in the View/Edit screen if AutoCert was triggered.
Performance	35571	Resolved Performance issue which frequently impacted the IVT workflow.
Scheduling	35572	Resolved Scheduling issue which caused errors to occur when attempting to access scanned documents.

SORTED BY CATEGORY AND REDMINE

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

