



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.078

UPDATED JULY 18, 2024

TABLE OF CONTENTS

Summary	2
Release Announcement.....	2
New Features.....	2
Resolved Issues.....	2
Feature Details	3
Patient Portal.....	3
Feature #35226 – Display message when loading online forms in Patient Portal	3
Version Details	4
Code Stream.....	4

PUBLICATION HISTORY

Revision	Author	Description
July 18, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2024.078 introduces better handling of Digital Forms within the Patient Portal and resolves Insurance Management, System performance and Scanning related issues.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Patient Portal	35226	Display message when loading online forms in Patient Portal	This enhancement to Patient Portal provides updated handling and presentation of Digital Forms.

SORTED BY CATEGORY AND REDMINE

Refer to the FEATURE DETAILS section below for configuration and usage information.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Insurance Management	35293	Resolved Insurance Management issue which impeded the printing of insurance forms following modifications to insurance data.
Patient Folder	35307	Resolved Patient Folder issue where changing the Order Type on multiple orders returned an error.
Performance	35380	Resolved Performance issue which caused memory leak errors when retrieving attachments.
Scanning	35378	Resolved Scanning issue where opening documents resulted in previously scanned documents being removed from the Attachment list.
System	34446	Updates System performance for appointment search times within the Patient and Referring portals.

SORTED BY CATEGORY AND REDMINE

FEATURE DETAILS

Patient Portal

Feature #35226 – Display message when loading online forms in Patient Portal

Summary

This enhancement to Patient Portal provides updated handling and presentation of Digital Forms.

Background

Previously, digital forms were not being presented automatically when `ORFEvaluatePendingDFsBySqlJob` (which causes digital forms logic to be evaluated via a job) was enabled, requiring users to either refresh the screen or switch to another tab and back to force the digital forms to appear.

Feature Description

With this change, RIS will check the jobs that evaluate pending questions every 10 seconds while loading appointment details. If a job has not completed, a message (customizable via `PPOnlineLoadFormsText`) will be displayed to the user indicating that the system is currently evaluating online forms.

Once the job has been executed successfully, the online form will be displayed.

Configuration Instructions

No System Administrator actions are necessary to enable this feature.

RIS Client

Changes to ParagraphConfig Lookup Table

- Optionally, adjust the text of the new `PPOnlineLoadFormsText`.

Setting	Default	Purpose
PPOnlineLoadFormsText	Evaluating pending or completed online forms to review.	Message to display to the patient when evaluating online forms by checking SQL job. Added in #35226.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

