



eRAD RIS

# RELEASE ANNOUNCEMENT

Build 4.2023.102

UPDATED NOVEMBER 30, 2023

## TABLE OF CONTENTS

**Summary** ..... 2

    Release Announcement .....2

    Resolved Issues .....2

**Version Details** ..... 3

    Code Stream .....3

## PUBLICATION HISTORY

Revision	Author	Description
November 30, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>▪ Commercial release.</li> </ul>

# SUMMARY

## Release Announcement

*This release of eRAD RIS 4.2023.102 resolves a collection of Digital Forms and other issues, as well as two service-only issues.*

## Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Database	32489	<b>This is a service-only change with no customer-facing functionality.</b> Enhanced Database performance by adding two missing indexes.
Digital Forms	34616	Resolved Digital Forms issue where the entire Date Type toolbar was erroneously hidden in the Digital Forms Editor.
Digital Forms	34702	Resolved Digital Forms issues where forms failed to load correctly if the Patient Alert code included apostrophe character.
Digital Forms – Drawing Tools	34576	<b>This is a service-only change with no customer facing functionality.</b> Addressed a Digital Forms Drawing Tool issue by reverting to the previous dll version.
Inbound Document	34648	Resolved Inbound Document problem where the Patient Search field was disabled for SecurePic.
Insurance Management	34659	Addressed Insurance Management problem where insurance card icons appeared blank following a reschedule.
Worklists	34200	Resolved Worklists problem where the All Signed By Date worklist was not populating correctly.

SORTED BY CATEGORY AND REDMINE

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:

