



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.077

UPDATED JULY 17, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement..... 2

 Resolved Issues..... 2

 New Known Limitations..... 2

Version Details 3

 Code Stream..... 3

PUBLICATION HISTORY

Revision	Author	Description
July 17, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2024.077 updates the CURES import tool to support CMS125v11 and resolves a C-CDA exporting issue.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
CCDA	35698	Resolved CCDA issue which prevented the C-CDA from being exported if an invalid NPI was present within the data.
CURES	35297	Updates CURES import QRDA tool to support CMS125v11.

SORTED BY CATEGORY AND REDMINE

New Known Limitations

The following new Known Limitations were identified with this release:

- **BUG #35640 – C-CDA PREVIEW AND EXPORT HTML FUNCTIONALITY IS NOT WORKING**
 - **Issue:** C-CDA preview, and export HTML functionality is hindered when PPCCDAReportHoldBackMessge system config contains line breaks.
 - **Impact:** An "unexpected error" displays when clicking [Preview](#). Users are unable to preview and export C-CDA documents in HTML format.
 - **Workaround:** Remove line breaks (
 tags) from the system config value.
- **FEATURE #35746 – C-CDA GENERATION WHEN MISSING NPI/SIGNING RAD (from #35698)**
 - In a future release, FEATURE #35746 – C-CDA GENERATION WHEN MISSING NPI/SIGNING RAD is planned to introduce a new workflow to manage any export failures logged due to the changes introduced with #35698.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

