



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.013

UPDATED FEBRUARY 27, 2025

TABLE OF CONTENTS

Summary	2
Release Announcement	2
Resolved Issues	2
Version Details	3
Code Stream	3

PUBLICATION HISTORY

Revision	Author	Description
February 27, 2025	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of *ERAD RIS 4.2025.013* resolves a scanned documents interface issue.



*This release resolves the Known Limitation identified in the **ERAD RIS SERVICE BULLETIN 20250219 - SCANNED DOCUMENTS MISMATCH**.*

This medium-severity issue was identified by the Development Team. In unusual circumstances, scanned documents may be incorrectly assigned. A software fix and data fix strategy are available.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Interfaces	36593	<p>Resolved Interfaces issue that caused data corruption by reassigning scanned documents incorrectly.</p> <p>This medium-severity issue was identified by the Development Team. In unusual circumstances, scanned documents may be incorrectly assigned.</p> <p>Data Fix</p> <p>A data fix strategy for previously impacted records is available. The Development Team has developed several database queries to identify affected records and will work with affected sites to prepare a combination of automated and manual updates to restore the majority of impacted records.</p>

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

