



eRAD RIS

RELEASE ANNOUNCEMENT

Rapid Release 4.2024.041.0.2

UPDATED APRIL 11, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 Resolved Issues2

Version Details 3

 Code Stream3

PUBLICATION HISTORY

Revision	Author	Description
April 11, 2024	Kevin Brooks / Hilary Saltmarsh	<ul style="list-style-type: none"> Commercial release.

SUMMARY

Release Announcement

This rapid release of ERAD RIS 4.2024.041.0.2 updates Provider Portal to resolve an issue with Single Sign-On (SSO) with Epic EMR systems.

Resolved Issues


This release resolves the following issues:

Category	Redmine #	Subject
Provider Security	35404	Resolved a Provider Portal issue where IVIS was not loading images from SSO generated encrypted links.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:



WARNING: *The functionality of this release has not yet been merged forward into any other release.*

Please carefully review the code stream diagram of any future releases before deploying to ensure that the functionality of this release will not be lost when upgrading.