



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.025

UPDATED MAY 2, 2023

PUBLICATION HISTORY

Revision	Author	Description
May 2, 2023	Michelle Mahabir, Kevin Brooks / Hilary Saltmarsh	<ul style="list-style-type: none">Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2023.025 resolves a collection of fixes to resolve errors and UI issues, including a Provider Portal issue which caused buttons to intermittently disappear from the Appointment screen.

New Features

No new features were introduced with this release.

Resolved Issues

This release resolves the following issues:

SORTED BY CATEGORY AND REDMINE

Category	Redmine #	Subject
Electronic Medical Records (EMR)	32598	Resolved and EMR issue where Document Distribution jobs for EMR orders went to Error status when there was not active referring address. Document Distribution will now create an EMR job in Pending status and a Print job in Error status.
Insurance Eligibility	8486	This is a service-only change with no customer-facing functionality.
Insurance Eligibility	33574	Resolved an Insurance Eligibility issue where the PreApproved Payment amount was not displayed in the ATC summary when the ATC was zero with no Prior Balance.
Localization	33395	Updated the UI message that appears to users attempting to change a patient alert when they do not have permissions to "You do not have permission to update patient alerts."
PACS	33432	Resolved PACS issue where Radiologists would encounter an error and RIS would freeze after signing a study from the Suspended WL.
Provider General Display	33671	Resolves a Provider Portal issue that was causing the Edit Order and Cancel Order buttons to intermittently disappear due to inconsistencies reading user configuration data.
Scanning	33630	Resolved a Scanning issue where adding attachments to the highest accession number in the Exams tab would return an error.
Worklists	32179	Resolved Worklists UI display issue where an extra tooltip popup with the hotkeys for every flag in the row would appear when hovering over flags.
Worklists	33389	Resolved Worklists issue where two columns for Last Contact Type were appearing on the Orders To Schedule Worklist.

New Known Limitations

No new known limitations were identified with this release.