



eRAD RIS

HOTFIX RELEASE ANNOUNCEMENT

Hotfix 4.2024.062.3

UPDATED JUNE 17, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement..... 2

 Resolved Issues..... 2

Version Details 3

 Code Stream..... 3

PUBLICATION HISTORY

Revision	Author	Description
June 17, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> ▪ Commercial release.

SUMMARY

Release Announcement

This hotfix of ERAD RIS 4.2024.062.3 resolves two issues pertaining to the generation of the C-CDA.

Resolved Issues

This release resolves the following issues:


Category	Redmine #	Subject
CCDA	35662	Resolved C-CDA issue which caused an error when generating a C-CDA due to a missing null check.
CCDA	35669	Resolved C-CDA issue where the active NPI was not being considered post person merge causing an error when generating a C-CDA.

SORTED BY CATEGORY AND REDMINE

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:



	Latest commercialized release
	Story release
	Hotfix release
	Rapid release
	Internal test release
	Recalled release