



eRAD RIS

RELEASE ANNOUNCEMENT

Rapid Release 4.2024.056.0.1

UPDATED MAY 8, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement..... 2

 New Features..... 2

Feature Details 3

 Patient Forms..... 3

 Service-Only Feature #35212 Support hosting Online Registration as a nested IIS application..... 3

Version Details 4

 Code Stream..... 4

PUBLICATION HISTORY

Revision	Author	Description
May 8, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> Commercial release.

SUMMARY

Release Announcement

This rapid release of ERAD RIS 4.2024.056.0.1 adds a service-only feature with new hosting options for the Online Registration application but has no customer-facing changes.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Patient Forms	35212	Support hosting online Registration as a nested IIS application	This is a service-only change with no customer-facing functionality.

Refer to the FEATURE DETAILS section below for configuration and usage information.

FEATURE DETAILS

Patient Forms

Service-Only Feature #35212 Support hosting Online Registration as a nested IIS application

Summary

This enhancement to Patient Forms adds support for hosting Online Registration as a nested IIS application on the Patient Portal server.

This is a service-only change with no customer-facing functionality.

Background

Previously, hosting Online Registration as a nested IIS application was not supported for the Patient Portal server. Selecting the `GO TO INCOMPLETE FORMS` button the Patient Portal would incorrectly generate the path at the parent level (patient portal), leading to the error message "Issue while loading digital form" in environments with the OR functioning as a nested IIS application.

Feature Description

With this change, a new `ORPortalCompleteURL` App Configuration setting allows specifying the full URL path to OnlineForms when the OR is deployed as a nested IIS application.

Service Team Upgrade Instructions

The Service Team must complete the following actions to deploy this feature:

Changes to applicationsettings.config

1. Configure the value for the new `ORPortalCompleteURL` (default is blank):
 - a. When Online Registration is NOT deployed as a nested IIS application, the value should be left blank, and the existing configuration value for `ORPEXternalBaseURL` will be used.
 - b. Otherwise, the value should specify the complete url:

`www.patientportal.com/{OnlineForms}/api/GetSSOURL`



For more detailed guidance on configuring digital forms within the portal, refer to the documentation for Feature #34140.


Configuration Instructions

Service Team assistance is required to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:



WARNING: *The functionality of this release has not yet been merged forward into any other release.*

Please carefully review the code stream diagram of any future releases before deploying to ensure that the functionality of this release will not be lost when upgrading.