



eRAD RIS

# CUSTOMER RELEASE NOTES

Hotfix Build 3.2022.8.29.6

UPDATED MARCH 24, 2023

# Table of Contents

<b>Summary .....</b>	<b>3</b>
Intended Audience .....	3
Release Content .....	3
Who Is Affected .....	3
<b>New Settings .....</b>	<b>4</b>
<b>New Features .....</b>	<b>5</b>
<b>Resolved Issues and Known Limitations .....</b>	<b>6</b>
Resolved Issues.....	6
New Known Limitations.....	6
<b>Version Details .....</b>	<b>7</b>
Package Contents .....	7
Code Stream .....	8
eRAD RIS Release Version Numbers.....	9

## Publication History

Revision	Author	Description
March 24, 2023	Kevin Brooks / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>▪ Commercial release.</li> </ul>

# SUMMARY

## Intended Audience

The intended audience for this CUSTOMER RELEASE NOTES document is the RIS Administration team for all eRAD RIS customers.

The intent of this document is to describe the content of the build or hotfix with sufficient detail for a Customer to be able to understand the features, and for the Customer's RIS Administrator to be able to perform routine configuration of features.

This information should be used by all customers to determine how these changes affect their organization's workflow, and to plan their organization's upgrade strategy for eRAD RIS.



*There is no separate SERVICE RELEASE NOTES edition of this document.*

## Release Content

*This hotfix release of eRAD RIS 3.2022.8.29.6 resolves an urgent Radiology Reporting issue displaying attachments when navigating to the next patient exam (for non-Powerscribe 360 users).*

## Who Is Affected

This version is recommended to be applied on eRAD RIS v3.2022.8.29.X installations.

# NEW SETTINGS

No settings were added or updated with this release.

# NEW FEATURES

No new features were introduced with this release.

# RESOLVED ISSUES AND KNOWN LIMITATIONS

## Resolved Issues

This release resolves the following issues:

Redmine #	Subject
33714	Resolved Radiology Reporting issue where displayed attachments were not refreshing when navigating to the next patient exam (for non-Powerscribe 360 users).


















## New Known Limitations

No new known limitations were identified with this release.

# VERSION DETAILS


## Package Contents

The release package includes the following folders:

 @Hotfixes	4/20/2012 8:51 AM
 _Documentation	3/22/2023 4:13 PM
 _ReleaseNotes	3/24/2023 3:31 PM
 Client Application	3/22/2023 4:13 PM
 DB	3/22/2023 4:13 PM
 External WebAPI	3/22/2023 4:13 PM
 Identity Service	3/22/2023 4:13 PM
 Management Reports	3/22/2023 4:13 PM
 PACS Citrix Bridge	3/22/2023 4:13 PM
 RIS Service	3/22/2023 4:13 PM
 Service Tools	3/22/2023 4:13 PM
 Web Digital Forms	3/22/2023 4:13 PM
 Web Patient Connect	3/22/2023 4:13 PM
 Web Referring Connect	3/22/2023 4:13 PM
 Web UM Connect	3/22/2023 4:13 PM
 Build_2022.8.29.6.zip	3/22/2023 4:13 PM
 RISServerMasterCert.pfx	3/31/2016 1:38 PM

## Code Stream

The following source code branches have been merged into this release:



LEGEND:

Light Green = Previously Released software

Gray = Internal version, non-release version

Bright Green = Current Release



**WARNING:** The functionality of this release has not yet been merged forward into any other release.

Please carefully review the code stream diagram of any future releases before deploying to ensure that the functionality of this release will not be lost when upgrading.



## eRAD RIS Release Version Numbers

The following table details the version identifiers for components in this release:

Build	Patch	UI Version	Core Version	WS Version	DB Version	Digital Forms	Patient Portal	UM Portal	Provider Portal	Notes
2022.8.1	-	3.22.8.1(3GB)	3.22.8.1	3.22.8.1	3.22.8.1.003455261	3.22.8.1	3.22.8.1	3.22.8.1	3.22.8.1	GUI, Web Services, DB, Portals and Digital Forms
2022.8.1	1	3.22.8.1.1(3GB)	3.22.8.1.1	3.22.8.1.1						GUI, Web Services,
2022.8.1	2	3.22.8.1.2(3GB)	3.22.8.1.2	3.22.8.1.2	3.22.8.1.203553632					GUI, Web Services and DB
2022.8.1	3			3.22.8.1.3			3.22.8.1.3	3.22.8.1.3	3.22.8.1.3	Portals only
2022.8.1	4	3.22.8.1.4(3GB)	3.22.8.1.4	3.22.8.1.4	3.22.8.1.403635561	3.22.8.1.4	3.22.8.1.4	3.22.8.1.4	3.22.8.1.4	GUI, Web Services, DB, Portals and Digital Forms
2022.8.29	-	3.22.8.29(3GB)	3.22.8.29	3.22.8.29	3.22.8.29.003564808	3.22.8.29.0	3.22.8.29.0	3.22.8.29.0	3.22.8.29.0	GUI, Web Services, DB, Portals and Digital Forms
2022.8.29	1	3.22.8.29.1(3GB)	3.22.8.29.1	3.22.8.29.1						GUI, Web Services,
2022.8.29	2	3.22.8.29.2(3GB)	3.22.8.29.2	3.22.8.29.2	3.22.8.29.203607010	3.22.8.29.2	3.22.8.29.2	3.22.8.29.2	3.22.8.29.2	GUI, Web Services, DB, Portals and Digital Forms
2022.8.29	3			3.22.8.29.2			3.22.8.29.3	3.22.8.29.3	3.22.8.29.3	GUI only, Web Serviced and Portals
2022.8.29	4	3.22.8.29.4(3GB)	3.22.8.29.4	3.22.8.29.4	3.22.8.29.403616642		3.22.8.29.4	3.22.8.29.4	3.22.8.29.4	GUI, Web Services, DB and Portals
2022.8.29	5	3.22.8.29.5(3GB)	3.22.8.29.5	3.22.8.29.5	3.22.8.29.503635705	3.22.8.29.5	3.22.8.29.5	3.22.8.29.5	3.22.8.29.5	GUI, Web Services, DB and Portals
2022.8.29	6	3.22.8.29.6(3GB)	3.22.8.29.6	3.22.8.29.6		3.22.8.29.6	3.22.8.29.6	3.22.8.29.6	3.22.8.29.6	GUI, Web Services, DB and Portals
2022.10.24	-	4.22.10.24(3GB)	4.22.10.24	4.22.10.24	4.22.10.24.003648552	4.22.10.24	4.22.10.24	4.22.10.24	4.22.10.24	GUI, Web Services, DB, Portals and Digital Forms
2022.10.24	1	4.22.10.24.1(3GB)	4.22.10.24.1	4.22.10.24.1	4.22.10.24.003648552	4.22.10.24.1	4.22.10.24.1	4.22.10.24.1	4.22.10.24.1	GUI, Web Services, DB, Portals and Digital Forms
2022.10.24	2	4.22.10.24.2(3GB)	4.22.10.24.2	4.22.10.24.2	4.22.10.24.203700209	4.22.10.24.2	4.22.10.24.2	4.22.10.24.2	4.22.10.24.2	GUI, Web Services, DB, Portals and Digital Forms
2022.10.24	3	4.22.10.24.3(3GB)	4.22.10.24.3							GUI only
2022.10.24	4	4.22.10.24.4(3GB)	4.22.10.24.4	4.22.10.24.4		4.22.10.24.4	4.22.10.24.4	4.22.10.24.4	4.22.10.24.4	GUI only, Web Serviced and Portals
2022.11.1	-	4.22.11.1(3GB)	4.22.11.1	4.22.11.1	4.22.11.1.003659860	4.22.11.1	4.22.11.1	4.22.11.1	4.22.11.1	GUI, Web Services, DB, Portals and Digital Forms
2022.11.2	-	4.22.11.2(3GB)	4.22.11.2	4.22.11.2	4.22.11.2.003709059	4.22.11.2	4.22.11.2	4.22.11.2	4.22.11.2	GUI, Web Services, DB, Portals and Digital Forms
2023.1	-	4.23.1(3GB)	4.23.1	4.23.1		4.23.1	4.23.1	4.23.1	4.23.1	GUI, Web Services, Portals and Digital Forms
2023.2	-	4.23.2(3GB)	4.23.2	4.23.2		4.23.2	4.23.2	4.23.2	4.23.2	GUI, Web Services, Portals and Digital Forms
2023.3	-	4.23.3(3GB)	4.23.3	4.23.3		4.23.3	4.23.3	4.23.3	4.23.3	GUI, Web Services, Portals and Digital Forms
2023.4	-	4.23.4(3GB)	4.23.4	4.23.4		4.23.4	4.23.4	4.23.4	4.23.4	GUI, Web Services, Portals and Digital Forms
2023.5	-	4.23.5(3GB)	4.23.5	4.23.5		4.23.5	4.23.5	4.23.5	4.23.5	GUI, Web Services, Portals and Digital Forms
2023.6	-	4.23.6(3GB)	4.23.6	4.23.6	4.23.006.0.003757897	4.23.6	4.23.6	4.23.6	4.23.6	GUI, Web Services, DB, Portals and Digital Forms
2023.7	-	4.23.7(3GB)	4.23.7	4.23.7	4.23.007.0.003780786	4.23.7	4.23.7	4.23.7	4.23.7	GUI, Web Services, DB, Portals and Digital Forms
2023.8	-	4.23.8(3GB)	4.23.8	4.23.8	4.23.008.0.003781109	4.23.8	4.23.8	4.23.8	4.23.8	GUI, Web Services, DB, Portals and Digital Forms
2023.9	-	4.23.9(3GB)	4.23.9	4.23.9	4.23.009.0.003786508	4.23.9	4.23.9	4.23.9	4.23.9	GUI, Web Services, DB, Portals and Digital Forms
2023.10	-	4.23.10(3GB)	4.23.10	4.23.10	4.23.009.0.003786508	4.23.10	4.23.10	4.23.10	4.23.10	GUI, Web Services, DB, Portals and Digital Forms
2023.11	-	4.23.11(3GB)	4.23.11	4.23.11	4.23.011.0.003791263	4.23.11	4.23.11	4.23.11	4.23.11	GUI, Web Services, DB, Portals and Digital Forms
2023.12	-	4.23.12(3GB)	4.23.12	4.23.12		4.23.12	4.23.12	4.23.12	4.23.12	GUI, Web Services, DB, Portals and Digital Forms