



eRAD RIS

RELEASE ANNOUNCEMENT

Rapid Release 4.2023.074.0.1

UPDATED SEPTEMBER 26, 2023

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 Resolved Issues2

Version Details 3

 Code Stream3

PUBLICATION HISTORY

Revision	Author	Description
September 26, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> Commercial release.

SUMMARY

Release Announcement

This rapid release of eRAD RIS 4.2023.074.0.1 resolves issues affecting the Provider Exam Detail Page usability.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Provider Exam Detail Page	34366	Addressed Provider Exam Detail Page issue which caused the report download button to become disabled if the exam page loaded slowly.
Provider Exam Detail Page	20645	Resolved Provider Exam Detail Page where bypassing the "I agree" validation prior to downloading a report led users to a blank webpage.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

