

# **Server Update**

**For eRAD RIS  
Version 1.0  
Build 38  
Update 38.007**

## **Table of Contents**

1. Purpose.....	3
Defects resolved.....	3
2. Package Contents .....	4
3. Intended Audience .....	4
4. Who is affected .....	4
5. eRAD RIS Release Version Numbers .....	5
6. Installing .....	6

## **RADNET, Inc. – eRAD RIS Server Update**

### **1. Purpose**

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update includes the following:

#### **Defects resolved**

Redmine #	Subject	Priority
2884	Hold Code Completed Status	Normal
2898	Verify and continue does not work	High
2904	HL7 messages are changing study status back in the workflow post-ExamDone	Urgent
2905	Reporting - Interactive dictation getting lost in data when refocusing after tab switch	High
2908	Reporting - inbound reporting may not open accession numbers tried once with an error	High

## **2. Package Contents**

Files included for this server update:

1. GUI ZIP Updates folder contains
  - a. rRIS\_1.2.38.15768.zip
2. RIS Service Updates folder contains
  - a. CDSJobManager.pdb
  - b. CDSJobManager.dll
  - c. CommonLib.pdb
  - d. CommonLib.dll
  - e. DDEngine.pdb
  - f. DDEngine.dll
  - g. rRISServices.pdb
  - h. rRISServices.exe
  - i. rRISServices.vshost.exe
  - j. ScheduledTaskManager.pdb
  - k. ScheduledTaskManager.dll
  - l. SignatureApprovalService.pdb
  - m. SignatureApprovalService.exe
  - n. WCFService.pdb
  - o. WCFService.dll
  - p. WebApi.pdb
  - q. WebApi.dll
  - r. Wedge.pdb
  - s. Wedge.exe

## **3. Intended Audience**

This document is created by the RIS Development team for the RadNet RIS Implementation team.

## **4. Who is affected**

Previous versions of build 38.

Server updates may not be cumulative; therefore this server update must be applied directly to the previous server update.

## 5. eRAD RIS Release Version Numbers

Build	Patch	UI Version	Core Version	WS Version	DB Version	Notes
37	-	1.2.37.12843	1.2.37.12782	1.2.37.12843	1.2.37.12843	
37	1	1.2.37.13040	1.2.37.13040	1.2.37.13039	1.2.37.13003	
37	2	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	
37	3	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	Just WebAPI updated, no version changes
37	4	1.2.37.13171	1.2.37.13040	1.2.37.13165	1.2.37.13171	
37	5	1.2.37.13347	1.2.37.13040	1.2.37.12979	1.2.37.12979	Bug was discovered with version numbers, this is why it appears we went backwards in the version number but we really did not.
37	6	1.2.37.13665	1.2.37.13040	1.2.37.13580	1.2.37.12979	
37	7	1.2.37.13807	1.2.37.13040	1.2.37.13799	1.2.37.13807	
37	8	1.2.37.13807	1.2.37.13040	1.2.37.13827	1.2.37.13807	
37	9	1.2.37.13807	1.2.37.13040	1.2.37.13839	1.2.37.13807	
37	10	1.2.37.13895	1.2.37.13040	1.2.37.13895	1.2.37.13895	
37	11	1.2.37.14161	1.2.37.13040	1.2.37.14161	1.2.37.14161	
37	12	1.2.37.14267	1.2.37.13040	1.2.37.14161	1.2.37.14267	
37	13	1.2.37.14326	1.2.37.13040	1.2.37.14326	1.2.37.14326	
38	-	1.2.38.14486	1.2.38.14486	1.2.38.14486	1.2.38.14486	
38	1	1.2.38.14748	1.2.38.14748	1.2.38.14748	1.2.38.14748	
38	2	1.2.38.15182	1.2.38.14748	1.2.38.14748	1.2.38.14748	
38	3	1.2.38.15251	1.2.38.14748	1.2.38.14748	1.2.38.14748	
38	4	1.2.38.15327	1.2.38.14748	1.2.38.15290	1.2.38.14748	
38	5	1.2.38.15513	1.2.38.14748	1.2.38.15513	1.2.38.15451	
38	6	1.2.38.15618	1.2.38.14748	1.2.38.15513	1.2.38.15618	
38	7	1.2.38.15768	1.2.38.14748	1.2.38.15768	1.2.38.15618	

## 6. Installing

- 1) Copy the zip file to be **rRIS\_1.2.38.15768.zip** file from the GUI ZIP Updates folder to the RISDeployment\Build folder. Move the rRIS\_1.2.38.15618.zip folder to its backup location. No changes to the ris.exe.config settings file that is contained within the .zip folder.
- 2) Stop each service listed below. Copy the files from the RIS Service Updates folder to the directory for each rRISService.  
*rRISService directories include:*
  - a) *rRISService*
  - b) *CDS rRISService*
  - c) *Wedge rRISService*
  - d) *DDE rRISService*

Restart each service.