



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.107

UPDATED OCTOBER 24, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 Resolved Issues2

 New Known Limitations2

Version Details 3

 Code Stream3

PUBLICATION HISTORY

Revision	Author	Description
October 24, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2024.107 resolves a recurring error dialog that appears within the Scheduling workflow.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Scheduling	35969	Resolved Scheduling issue which caused a recurring error dialog throughout the scheduling workflow.

New Known Limitations

The following significant limitations have been identified and should be communicated to affected users:

- **BUG #36114 PERFORMANCE ISSUES USING REPORT TOGETHER PROCEDURE PLANS (FROM #35969)**
 - **Issue:** Procedure plans that have the report together flag (which assigns a linked id to procedures in the plan) are taking a long time to complete transactions on systems with high volumes of records.
 - **Impact:** Response times of several minutes have been observed when performing scheduling, rescheduling, arriving, perform exam and other actions on report together procedures. Systems with lower record volumes may experience little to no issues.
 - **Workaround:** No known workaround. A fix is planned for a future release. Note that this is a pre-existing issue that was also observed in prior releases.



Carefully review all Known Limitations identified here and in the release notes of all other versions being applied and communicate relevant items to all those affected.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

