



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.046

UPDATED SEPTEMBER 5, 2025

TABLE OF CONTENTS

Summary	2
Release Announcement.....	2
New Features.....	2
Resolved Issues.....	2
Feature Details	3
Feature #35739 – Support Digital Forms Required Fields Attribute for Guided Workflow.....	3
Version Details	4
Code Stream.....	4

PUBLICATION HISTORY

Revision	Author	Description
September 5, 2025	Michelle Mahabir/ Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2025.046 enhances Guided Workflow support of required fields on Digital Forms and resolves various Guided Workflow scheduling issues.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
GW-Mammo	35739	Support Digital Forms Required Fields Attribute for Guided Workflow	This enhancement to Guided Workflow - Mammo adds support for required fields on Digital Forms by introducing an attribute to properly enforce required questions during the workflow.

SORTED BY CATEGORY AND REDMINE

Refer to the FEATURE DETAILS section below for configuration and usage information.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
GW - Mammo	37528	Resolved GW - Mammo issue which presented silent errors when accessing the Perform Exam screen after logging in.
GW - Mammo	37534	Resolved GW - Mammo issue that prevented progression through the workflow due to the Historical Carrier Selection process.
GW - Mammo	37553	Resolved GW - Mammo issue where Site-Level Procedure Plans were not being recognized within the Guided Workflow.

SORTED BY CATEGORY AND REDMINE

FEATURE DETAILS

Feature #35739 – Support Digital Forms Required Fields Attribute for Guided Workflow

Summary


This enhancement to Guided Workflow for Mammography adds support for required fields on Digital Forms by introducing an attribute to properly enforce required questions during the workflow.

Background

Previously, in Guided Workflow, digital form validation did not account for UI actions or workflow status. This caused all form questions to be required immediately, rather than only at the intended step (e.g., Check-In), which led to questions appearing too early in the process.

Feature Description

The Digital Form Editor now includes a `Guided Workflow Required` checkbox. Checking this box ensures the form will be required during Guided Workflow.



The screenshot displays the Digital Form Editor interface. At the top, there are several checkboxes: 'Save answers at the order level' (checked), 'Copy forms completed with the past 0 days' (unchecked), 'Prefill answers entered within the past 0 days' (unchecked), 'Right to left' (unchecked), 'Delay expression evaluation' (unchecked), and 'Delete hidden answers' (unchecked). Below these, there is a section for 'Body' with the text 'Are you experiencing any NE'. The 'Answer Type' is set to 'Radio Button'. The 'Short name' is 'any_new_issues'. There is a 'Map question' button and a 'Map' button. Below this, there are checkboxes for 'Patient Required' (unchecked), 'RIS Required' (unchecked), 'On UI Action' (unchecked), 'Edit' (unchecked), 'On or After Status' (unchecked), and 'Guided Workflow Required' (checked). At the bottom, there is a 'Font size' dropdown, a 'Bold' button, a 'Font color' dropdown set to 'Default', and a 'Background color' dropdown set to 'Default'. A red arrow points to the 'Guided Workflow Required' checkbox.

Configuration Instructions

No System Administrator actions are necessary to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

