



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.058

UPDATED NOVEMBER 19, 2025

TABLE OF CONTENTS

Summary	2
Release Announcement	2
New Features	2
Resolved Issues	2
Feature Details	3
Digital Forms	3
Feature #37602 - Limit SmartOnboarding Digital Form access to active date studies and hide header/footer elements	3
Version Details	4
Code Stream	4

PUBLICATION HISTORY

Revision	Author	Description
November 19, 2025	Kevin Brooks, Michelle Mahabir	▪ Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2025.058 limits SmartOnboarding Digital Form access to active date studies and resolves a RADAR connectivity issue to the SignalR service.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject
Digital Forms	37602	Limit SmartOnboarding Digital Form access to active date studies and hide header/footer elements.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
RADAR	36323	Resolved RADAR issue that disrupted connectivity to the SignalR service.

FEATURE DETAILS

Digital Forms

Feature #37602 - Limit SmartOnboarding Digital Form access to active date studies and hide header/footer elements

Summary

This enhancement to Digital Forms updates questionnaires to only display when a study is scheduled for the current day, as required for the Smart Onboarding workflow.

Additionally, the form's header, footer, and `Submit` button can now be suppressed when the questions are displayed on a smart onboarding device.

Configuration Instructions

No System Administrator actions are necessary to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

