



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.061

UPDATED DECEMBER 3, 2025

TABLE OF CONTENTS

Summary	2
Release Announcement	2
New Features	2
Feature Details	3
Insurance Management	3
Feature #37881 - Config to Disable Cascade Removal of Expired Carriers Functionality	3
Version Details	4
Code Stream	4

PUBLICATION HISTORY

Revision	Author	Description
December 3, 2025	Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2025.061 introduces a new system configuration option to disable the "Cascade Removal of Expired Carriers to Other Visits" functionality introduced in v4.2025.057.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Insurance Managment	37881	Disable Cascade Removal of Expired Carriers Functionality	This enhancement to Insurance Management introduces a new system configuration option to disable the "Cascade Removal of Expired Carriers to Other Visits" functionality introduced in V4.2025.057.

Refer to the FEATURE DETAILS section below for configuration and usage information.

FEATURE DETAILS

Insurance Management

Feature #37881 - Config to Disable Cascade Removal of Expired Carriers Functionality

Summary

This enhancement to Insurance Management introduces a new system configuration option to disable the "Cascade Removal of Expired Carriers to Other Visits" functionality introduced in v4.2025.057.

Background

Previously, FEATURE #36718 updated Insurance Management carrier handling in RIS to remove an expired carrier from all other visits when the carrier has been marked as expired via the RIS UI. However, it did not provide a configuration option to disable the new functionality. When used in production, users found that it was interfering with current worklist filtering causing agents to lose track of precert work already in process.

Feature Description

With this change, a new `AutomaticRemovalOfExpiredCarriersFromRelatedVisits` System Configuration setting allows this workflow to be disabled.

Known Limitations

The following significant limitations have been identified and should be communicated to affected users:

- **BUG #37888 - EXPIRED CARRIER EXAMS NOT GETTING PLACED ON THE IVT WORKLIST**
 - **Issue:** When FEATURE #36718 was introduced, to remove an expired carrier from other visits, it did not accommodate those items that did not initially require authorization but then needed insurance verification because the carrier was expired on a different exam in the folder.
 - **Impact:** Affected items are not going to the IVT WL without intervention.
 - **Workaround:** Manually re-opening and re-saving affected exams will move them to the IVT WL.

Configuration Instructions

System Administrators must complete the following actions to enable this feature:

RIS Client

Changes to SystemConfig Lookup Table

- Enable this feature by setting `AutomaticRemovalOfExpiredCarriersFromRelatedVisits` to `True` and configuring per FEATURE #36718 introduced in ERAD RIS CUSTOMER RELEASE NOTES 4.2025.057.

The following related settings were added or updated:

Setting	Default	Purpose
AutomaticRemovalOfExpiredCarriersFromRelatedVisits	Value=Boolean, Set to [False] to disable, Default=[False]	When [True], expired carriers are automatically removed from other visits when not yet scheduled or status is Scheduled (per Feature #36718). Updated by #37881

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

